
# PREMIUM BRACING AND PERFORMANCE APPAREL





# Table of Contents

Message from the CEO
Premium Materials
Stabilization5Knee5Range of Motion5Aluminum Hinge7Conforma™ Self Shaping Hinge9Spine11Wrist13
Foot & Ankle
Compression.       19         Knee.       21         Calf.       23         Spine.       24         Thigh & Hip       24         Wrist, Elbow, Foot & Ankle       25         Apparel       26
Quick Reference Guide

Dear Customer,

Bio Skin was born by fusing superior fabrics together to create a powerful material—thin, breathable, hypoallergenic, and compressive. More than 20 years later, Bio Skin offers an entire portfolio of products with innovative designs and the highest quality textiles.

Our success as a brand is the result of serving customers who demand innovation. Bio Skin offers premium braces to complement the in-office care that patients receive in premium clinics.

Bio Skin products embody the finer aspects of bracing including muscle expansion during flexion, sensitivity of the popliteal area behind the knee, and optimal compression needed for recovery and performance.

We are committed to preventing injury, accelerating recovery, and improving human performance. We continue to discover, develop, and deliver products that enhance quality of life.

Thank you for including us in your pursuit to Be Better™.

- Tanner Cropper, CEO






# Premium Materials

We've designed premium materials with body climate control that are hypo-allergenic, breathable, light, thin, and strong. These special materials form the foundation upon which we build our entire product line and allow us to provide your patients with the best and most comfortable braces on the market.

## **Ultima**<sup>1</sup>

#### Ultra-Thin, High Compression



Ultima is our original, tri-laminate material. Ultima is amazingly thin and provides high-level compression with powerful Lycra and a polyurethane film. The film allows moisture vapor transfer so perspiration evaporates, facilitating the body's natural temperature regulation system. The Ultima material is used widely in our product line.

## Stratus'

### Extra Breathability, Forgiving Fit



Stratus has more body than Ultima for people who like the feeling of extra support. Stratus is a laminated material with high-power Lycra on the outsides and a porous membrane in the middle. The more elastic membrane in Stratus allows it to provide more dynamic compression and a more forgiving fit.

### 2SL<sup>™</sup>

#### Front Closure Technology



2SL features a Velcro®-compatible surface to provide compression and flexibility in anchoring and closures. The 2SL material is used primarily in our Front-Closure Knee Skins to provide compression and allow total adjustability.

## Cirrus<sup>™</sup>

#### Single Layer, High Performance Stretch



Cirrus is made of a single layer of high-performance Lycra. Cirrus is incredibly thin and perfectly breathable, while providing a moderate level of compression.



# RANGE OF MOTION HINGED KNEES





Adjustable Range of Motion Hinge for Gladiator

## Gladiator<sup>™</sup> Knee Ligament Protection

The Gladiator uses a rigid hinge to protect knee ligaments, including post-operatively. The rigid strapping system firmly anchors the Gladiator to the leg, avoiding the need for larger, more expensive braces. Hinges are easily set to limit flexion and/or extension.

#### **FEATURES**

- Length: 15"
- Silicone gel patella ring (pull-on only)
- Easily adjustable flexion / extension settings without tools
- Open patella
- Pull-on or front-closure design
- · Fits left or right

#### INDICATIONS

- ACL
- PCL
- LCL
- MCL
- Knee Instability

Product	Product #
Gladiator - Open Patella - Ultima™ Material	4134*
Gladiator - Open Patella - Stratus™ Material	6136*
Gladiator - Front Closure - Open Patella	5136*


**Note:** To indicate size, **replace \*** with corresponding number from the sizing chart on the back cover.

# CrossFire<sup>™</sup> Line

CrossFire products offer two major differences compared to other products. First, all CrossFire products have a unique back panel design that offers comfortable compression throughout the full range of knee flexion. Second, each CrossFire product is available in a front closure option we call Hybrid<sup>®</sup>, because it offers the compression of a pull on sleeve and the ease of donning and doffing of a wrap around.

## ActiveSeam<sup>™</sup> Stitching

With no overlapping material, our seams are flatter than other stitches and are 30% stronger than standard seams.



### Hybrid<sup>™</sup> Front Closure

An inner sleeve helps patients with donning and provides better compression than any wrap-around on the market.



#### Cirrus<sup>™</sup> Our Thinnest Material

The CrossFire cut uses our thinnest and most elastic material – offering optimal compression without bunching or binding.









Adjustable Range of Motion Hinge for Q Lok ROM

## O Lok ROM™ Patella Traction System Featuring the New SlimLok™ Adjustable Hinge

The Q Lok ROM is designed to diminish the effects and aid in the rehabilitation of patellofemoral pain. The Dynamic Patella Traction System exerts a medial pull on the patella to improve patellar tracking in the trochlear notch and increase surface contact area to mitigate pain.

#### **FEATURES**

- Length: 13"
- Addresses glide and tilt
- Fits left or right
- Easily adjustable flexion / extension settings
- ActiveSeam stitch

#### **INDICATIONS**

- Anterior Knee Pain
- · Maltracking Patella
- Patellofemoral Osteoarthritis
- Chondromalacia

Product	Product #
Q Lok ROM - Closed Patella	6166*
Q Lok ROM - Closed Patella - Hybrid™ Front Closure	5166*

Note: To indicate size, replace \* with corresponding number from the sizing chart on the back cover.



## CrossFire ROM™ Knee Ligament Protection Featuring the New SlimLok™ Adjustable Hinge

The CrossFire ROM is designed to protect knees with ligament tears or sprains. The rigid uprights are slightly flared to contour to the leg's natural shape and provide stability so patients can move and be active with confidence. The thin hinge is easy to adjust for flexion and extension control during patient rehabilitation.

#### **FEATURES**

- Length: 13"
- Fits left or right
- Easily adjustable flexion / extension settings
- ActiveSeam stitch

#### **INDICATIONS**

- Ligament Sprains / Tears
- Meniscus Tears
- Knee Instability

Product	Product #
CrossFire ROM - Closed Patella	6180*
CrossFire ROM - Closed Patella - Hybrid™ Front Closure	5180*

# Aircraft-grade Aluminum Hinges

Bio Skin combines superior compression textiles with aircraft grade aluminum in the Hinged Knee Skin line to bring you premium hinged knee support.

Both the CrossFire and Hinged Knee Skin are made with 6061 aluminum alloy. High strength-to-weight characteristics make this grade of aluminum the alloy of choice for airplane wings and the chassis of high-end sports cars. Bio Skin aluminum hinges are designed with smooth gears that match knee biomechanics with bi-centric rotation throughout flexion.



## CrossFire<sup>™</sup> with Aluminum Hinge

The CrossFire embodies our commitment to the importance of proper compression with a new design that maximizes style, comfort and stability. CrossFire knee braces use an ultra-thin textile and moderate compression to stabilize the knee and control swelling.

#### **FEATURES**

- Length: 13"
- Fits left or right
- ActiveSeam stitch
- Pull-on or front-closure design

#### **INDICATIONS**

- Ligament Sprains / Tears
- Meniscus Tears
- Knee Instability

Product	Product #
CrossFire Aluminum - Pull on w/ Standard Hinge - Closed Patella	6186*
CrossFire Aluminum - w/ Standard Hinge - Hybrid™ Front Closure - Closed Patella	5186*

**Note:** To indicate size, **replace \*** with corresponding number from the sizing chart on the back cover.



# Hinged Knee Skin<sup>™</sup> with Aluminum Hinge

This is the original Hinged Knee Skin—our most popular hinged knee brace. The Hinged Knee Skin is sleek, functional, durable, and very compressive.

#### **FEATURES**

- Offset hinge reduces pressure at the knee joint
- Strapping system minimizes migration
- Open patella
- Pull-on or front-closure design
- · Fits left or right

#### INDICATIONS

- LCL Tear or Sprain
- MCL Tear or Sprain
- Meniscus Tear

Product	Product #
Hinged Knee Skin Aluminum - Open Patella - Ultima™ Material	4130*
Hinged Knee Skin Aluminum - Open Patella - Stratus™ Material	6130*
Hinged Knee Skin Aluminum - Front Closure	5133*



# CONFORMA™ HINGED KNEES





## O Lok<sup>™</sup>

Shown to affect patellar positioning in weight-bearing MRI studies. Our unique design decreases pain by increasing surface contact area in the joint.

#### **FEATURES**

- · Addresses glide and tilt
- Closed patella
- Lightweight: 7 ounces
- Pull-on or front-closure design
- · Fits left or right

#### **INDICATIONS**

- Patellar Tilt and Glide
- Patellar Subluxation
- Post Surgical Lateral Release
- Patellofemoral Osteoarthritis

Product	Product #
Q Lok Dynamic Patella Traction Brace - Stratus™ Material	6164*
Q Lok Dynamic Patella Traction Brace - Front Closure	5164*

Note: To indicate size, replace \* with corresponding number from the sizing chart on the back cover.

# CrossFire<sup>™</sup> with Conforma<sup>™</sup> Hinge

Designed for active knees in need of extra support due to ligament tears or sprains. Each bicentric, nylon-polymer hinge weighs only one ounce, making this brace very light and comfortable.

#### **FEATURES**

- Length: 13"
- Fits left or right
- ActiveSeam stitch
- Pull-on or front-closure design

#### **INDICATIONS**

- Ligament Sprains / Tears
- Meniscus Tears
- Knee Instability

Product	Product #
CrossFire Pull on w/ Conforma Hinge - Closed Patella	6185*
CrossFire w/ Conforma Hinge - Hybrid™ Front Closure - Closed Patella	5185*

Note: To indicate size, replace \* with corresponding number from the sizing chart on the back cover.

# Self-shaping Conforma™ Hinge

Built with injection-molded nylon, the Conforma Hinge lives up to its name by conforming to the natural shape of the leg without any adjustment. This is the perfect solution for patients who want stability and comfort.



# **E** Lightweight

Each hinge weighs only one ounce, making the total weight for an entire brace a mere 6.75 ounces.



#### Low-Profile

The Conforma hinge is less than 5mm thick and the uprights are less than 2.5mm thick.



#### **Flexible**

Nylon plastic uprights conform perfectly to the shape of any leg, without the need for any alterations by the fitter.

# Visco Knee Skin<sup>™</sup> with Conforma<sup>™</sup> Hinge

Silicone ring reduces swelling around the joint.



- *Ideal after knee arthroscopy*
- Finger loops enable easy pull-on application
- Closed patella
- Fits left or right

#### **INDICATIONS**

- Knee Sprains and Strains
- Knee Joint Swelling
- · Meniscus Tear
- Osteoarthritis
- Patellar Tendonitis



Pull-on Front-closure

# Q Brace<sup>™</sup> with Conforma Hinge

Applies dynamic, lateral pressure to the patella.

#### **FEATURES**

- Open patella
- Pull-on or front-closure design
- Fits left or right

#### **INDICATIONS**

- Anterior Knee Pain
- Patellar Subluxation
- · Lateral Patella Tracking
- Medial/Superior Patella Tracking
- Medial Patella Tracking
- Medial/Inferior Patella Tracking

Product	Product #
Q Brace w/ Conforma Hinge	4136*
- Open Patella - Ultima™ Material	
Q Brace w/ Conforma Hinge - Front Closur	e 5138*

Note: To indicate size, replace \* with corresponding number from the sizing chart on the back cover.

#### Product # Visco Knee Skin with Conforma Hinge 6137 - Closed Patella - Stratus™ Material

Note: To indicate size, replace \* with corresponding number from the sizing chart on the back cover.

Pull-on Only

Pull-on Only

# Patella Stabilizer with Conforma Hinge

Stabilizes the patella and provides general support.

#### **FEATURES**

- Foam buttress around patella for stability
- Open patella
- Fits left or right

#### **INDICATIONS**

- Patella Stabilization
- · Knee Sprains and Strains
- Knee Joint Swelling
- Meniscus Tear
- Osteoarthritis
- Patellar Tendonitis

Product	Product #
Patella Stabilizer w/ Conforma Hinge	4132*
- Open Patella - Ultima™ Material	
Patella Stabilizer w/ Conforma Hinge	6132*
- Open Patella - Stratus™ Material	

Note: To indicate size, replace \* with corresponding number from the sizing chart on the back cover.

# Premium J<sup>™</sup> with Conforma Hinge

Elastic straps cross the knee to help seat the patella in the trochlear notch.

#### **FEATURES**

- Adjustable, lateral patellar tension straps
- Closed patella
- · Specify left or right

#### **INDICATIONS**

Pull-on Only

- Anterior Knee Pain
- · Lateral Tracking Patella
- Patellar Subluxation
- Patellofemoral Knee Pain

Product	Product #
Premium J with Conforma Hinge	6138*L
- Closed Patella - Stratus Material - Left	
Premium J with Conforma Hinge	6138*R
- Closed Patella - Stratus™ Material - Right	

Our Spine Line offers the world's best non-invasive solution for acute and chronic back pain. Only Bio Skin products offer Targeted Vector Compression to provide extra support to the posterior aspect of the spine.

#### **FEATURES**

- Targeted Vector Compression
- · Two hot/cold gel pads
- Easy-to-pull bilateral cinching system
- Newly designed cord cover



#### **INDICATIONS**

- Herniated Discs
- Degenerative Disc Disease
- Spondylolisthesis
- Osteoarthritis

# Baja<sup>™</sup> 627

#### **FEATURES**

 Posterior panel provides support from L-5 to L-1

#### **INDICATIONS**

• Lumbar Sprains and Strains

Product	Product #
Baja 627	3750*
Baja 627 - Trim to Fit - Universal Size	37504U

Note: To indicate size, replace \* with corresponding number from the sizing chart on the back cover.



# Baja<sup>™</sup> 631

#### **FEATURES**

 Posterior panel provides support from sacrum to T-9

#### **INDICATIONS**

• Radicular Pain

Product	Product #
Baja 631	3770*
Baja 631 - Trim to Fit - Universal Size	37704U

**Note:** To indicate size, **replace \*** with corresponding number from the sizing chart on the back cover.



Lumbar / Thoracic with Sagittal Control

# Baja<sup>™</sup> 637

#### **FEATURES**

- Posterior panel provides support from sacrum to T-9
- Removable side panels
- Sagittal control

#### **INDICATIONS**

- Radicular Pain
- · Lumbar Spinal Stenosis
- Post-Surgical Stabilization

Product	Product #
Baja 637	3760*
Baja 637 - Trim to Fit - Universal Size	37604U





## Back Flex™

Offers a flexible, posterior panel that supports the lumbar spine.

#### **FEATURES**

- Targeted Vector Compression
- Two hot/cold gel pads
- Easy-pull bilateral cinching system
- Flexible posterior panel

#### **INDICATIONS**

- Herniated Discs
- Degenerative Disc Disease
- Lumbar Sprains & Strains
- Spondylolisthesis
- Osteoarthritis

Product	Product #
Back Flex	3730*

**Note:** To indicate size, **replace \*** with corresponding number from the sizing chart on the back cover.



## Back Plus™

Offers a flexible, posterior panel that supports up to T-9 in the thoracic spine.

#### **FEATURES**

- Targeted Vector Compression
- Two hot/cold gel pads
- Easy-pull, bilateral cinching system Post-Kyphoplasty
- Shoulder straps
- Flexible posterior panel

#### **INDICATIONS**


- Compression Fractures
- Post-Vertebroplasty
- · Post-SCS Implant
- Post-IT Pump Implant
- Post-Fusion

Product Product # Back Plus

 $\textbf{Note:} \ \text{To indicate size,} \ \textbf{replace *} \ \text{with corresponding number from the sizing chart on the back cover.}$ 

# How Our Spine Line is Different

## Vector<sup>™</sup> Compression



Our Targeted Vector Compression focuses support directly on the posterior aspect of the spine. The result is extra pain relief and increased stability.

# Integrated Hot/Cold Therapy



Easy-to-use gel pads slide into the lumbar panel to offer extra pain relief.



DP2 - 6 1/2"

DP3 - 8 1/2"

# DP2<sup>™</sup> & DP3<sup>™</sup> Cock-Up Wrist

L3908

#### **FEATURES**

- Low profile and comfortable
- Allows full use of fingers and thumbs
- Removable/bendable aluminium dorsal and palmar stays
- Specify left or right

#### **INDICATIONS**

- Carpal Tunnel Syndrome
- Arthritis
- Wrist Sprain
- Bursitis

Product	Product #
DP2 Cock-Up Wrist Splint 6 1/2" - Left	5423*L
DP2 Cock-Up Wrist Splint 6 1/2" - Right	5423*R
DP3 Cock-Up Wrist Splint 8 1/2" - Left	5424*L
DP3 Cock-Up Wrist Splint 8 1/2" - Right	5424*R

Note: To indicate size, replace \* with corresponding number from the sizing chart on the back cover.





#### **FEATURES**

- Controls motion of the CMC and MP joints
- · Cock-up wrist support
- · Soft fleece material against skin
- · Adjustable aluminium stays
- Specify left or right

#### **INDICATIONS**

- · De Quervains
- Scaphoid Fracture
- Bursitis
- · Sprains and Strains

Product	Product #
Wrist/Thumb Spica - Left	5452*L
Wrist/Thumb Spica - Right	5452*R

Note: To indicate size, replace \* with corresponding number from the sizing chart on the back cover.



# Thumb Spica




#### **FEATURES**

- Controls motion of the CMC and MP joints
- Low profile
- · Soft fleece material against skin
- Fits left or right

#### **INDICATIONS**

- De Quervains
- Game Keeper's Thumb
- · Sprains and Strains

Product	Product #
Thumb Spica	5450*





TriLOK External Ligament System





## Subtalar Joint Control

- Valgus & Varus Correction
- Ankle Stabilization

### Lateral Column Unloading

- Ankle Stabilization
- Peroneal Tendon Dysfunction

### Medial Column Unloading

- Posterior Tibial Tendon Dysfunction
- Plantar Fasciitis



#### **FEATURES**

- Compressive base sleeve improves proprioception
- Easily fits into most shoes and sandals
- · Allows normal range of motion
- No laces

#### **INDICATIONS**

- Ankle Sprains
- Peroneal Tendinitis / Tendinosis
- · Chronic Ankle Instability
- Posterior Tibial Tendon Dysfunction
- Plantar Fasciitis / Heel Pain

Product Product # TriLok 5360\*

Note: To indicate size, replace \* with corresponding number from the sizing chart on the back cover.

Patent # 9,364,363 B2




# The Defining Difference

When ligaments and tendons of the foot are injured, the FootLok Strap takes over the responsibility to stabilize and support.



The white FootLok Strap wraps around the midfoot and forefoot and anchors to the medial or lateral side of the foot to control motion. The non-stretch FootLok Strap acts as an external ligament to provide suspension and support to weak and injured tendons and ligaments.







The AFTR<sup>™</sup> Line was named by combining the words "Ankle Foot Trauma Rehabilitation". It combines stabilization with high-level compression unique to Bio Skin®.

#### **FEATURES**

- Reduce Pain
- Improve Proprioception
- Increase Range of Motion
- Reduce Recovery Time







#### **INDICATIONS**

- Peroneal / Post-Tib Tendon Repair
- Lateral Ankle Repair
- Fusions: Midfoot, Hindfoot, Ankle
- Total Ankle Replacement
- Ankle Sprains
- Bunionectomy
- Forefoot Procedures
- Hindfoot Arthritis

## AFTR DC with Gel<sup>™</sup> (Dorsal Closure)

#### **FEATURES**

- Accommodates changes in foot & ankle volume
- Easy to apply during early, post-op phase
- Includes rigid stirrup strap
- Includes medial and lateral malleolar gel pads
- Trimmable in length and width 🧩
- Available with Figure 8 Compression Wrap (see below)

Product	Product #
AFTR DC with Gel	5386*
AFTR DC with Gel w/Figure 8 Compression Wrap (see below)	5386*CW

**Note:** To indicate size, **replace \*** with corresponding number from the sizing chart on the back cover.

## AFTR™ Pull-on

#### **FEATURES**

- Pull-on design
- •Trimmable in length ¾
- Includes rigid stirrup strap
- Available with Figure 8 Compression Wrap (see below)

Product	Product #
AFTR Pull-on	5380*
AFTR with Gel Pull-on (includes Malleolar Gel Pads)	5383*

**Note:** To indicate size, **replace \*** with corresponding number from the sizing chart on the back cover.

# AFTR w/Figure 8 Compression Wrap

#### **FEATURES**

- Figure 8 Compression Wrap
- Same features as AFTR Pull-on

Product	Product #
AFTR DC with Gel w/Figure 8 Compression Wrap	5386*CW
AFTR with Gel Pull-on w/Figure 8 Compression Wrap	5383*CW
AFTR Pull-on w/Figure 8 Compression Wrap	5380*CW



Anchor for: Hallux Control Strap and Weil Osteotomy Strap.

# Calibrated Midfoot Compression Wrap

29540

Provides an attachment point for the Weil Digit Strapping System without covering the ankle.

Made of hypo-allergenic, non-slip material to prevent migration.

Easily fits inside normal footwear.

Product	Product#
CMC Wrap 10-Pack -Universal Size	95362



Shown with Calibrated Midfoot Compression Wrap (sold separately). Also compatible with AFTR DC with Gel.

# HCS<sup>™</sup> (Hallux Control Strap<sup>™</sup>)

29550

Splinting - Maintains the hallux in a neutral position during critical, soft-tissue healing after bunionectomy.

**Stretching** - Stretches extensor tendons after primary healing to improve ROM and allow normal gait.

**Strengthening** - Provides resistance for patients to strengthen FHL for strong toe-off.

Note: Must be worn with AFTR DC with Gel, or the Calibrated Midfoot Compression Wrap.

Product	Product #
HCS 10-Pack	9534*T

**Note:** To indicate size, **replace \*** with corresponding number from the sizing chart on the back cover.



Shown with <u>Calibrated Midfoot Compression Wrap</u> (sold separately). Also compatible with <u>AFTR DC with Gel</u>, <u>AFTR with Gel</u> and <u>AFTR</u>.

# Weil Osteotomy Strap™

29550

Helps prevent dorsal drift (floating toe) after Weil Osteotomy by holding the toe in a neutral or plantarflexed position.


May be worn inside a walking boot, post-op shoe, or normal footwear and at night to plantarflex the digit and provide extensor stretch.

May be used non-operatively to treat metatarsalgia or plantar plate tears.


**Note:** Must be worn with any AFTR brace, or the Calibrated Midfoot Compression Wrap.

Product	Product #
Weil Osteotomy Strap 10-Pack -Universal Size	95372





# **COMPRESSION**



# Ten-7 The Knee Brace That DOESN'T Go On The Knee



#### **FEATURES**

- · Thigh Design
- Easily Hidden Under Clothing

#### INDICATIONS

- Chronic Knee Pain
- Patellofemoral Pain
- Osteoarthritis
- Meniscus Tear + Pain

Product	Product #
Ten-7 Pull-on	5146*

 $\textbf{Note:} \ \, \textbf{To indicate size, replace *} \ \, \textbf{with corresponding}$ number from the sizing chart on the back cover.



See product page for testimonials.



"I thought it was pain that I was going to have to live with for the rest of my life... Within 5 weeks, the pain was, for the most part, gone!"

Mike:: Tennis Pro



"I suffered a tear in the outside of my meniscus...within a couple of days of trying on the Ten-71 experienced a remarkable decrease in pain!"

Melissa:: Runner

Pull-on Front-closure

# Q Brace<sup>™</sup>

This unique, multi-functional, patella-tracking brace offers multidirectional tracking to the patella.

#### **FEATURES**

- Achieves greater dynamic force as the knee is flexed
- Open patella
- Pull-on or front-closure design
- Fits left or right
- · Available with a hinge for increased stability (see pg. 10)

#### **INDICATIONS**

- · Anterior Knee Pain
- Patellar Subluxation
- Lateral Patella Tracking
- Medial/Superior Patella Tracking
- · Medial Patella Tracking
- Medial/Inferior Patella Tracking

Product	Product #
Q Brace - Open Patella - Ultima™ Material	4161*
Q Brace - Front Closure	5163*

Note: To indicate size, replace \* with corresponding number from the sizing chart on the back cover.



Pull-on Only

Pull-on Only

## Visco Knee Skin<sup>™</sup>

Silicone ring reduces swelling and provides passive massage to control edema.

#### **FEATURES**

- Ideal after knee arthroscopy
- Pull-on finger loops
- Closed patella
- · Fits left or right
- Available with a hinge (see pg. 10)

#### **INDICATIONS**

- Knee Sprains and Strains
- Knee Joint Swelling
- Meniscus Tear
- Osteoarthritis
- Patellar Tendonitis

Product	Product #
Visco Knee Skin - Closed Patella	4170*
- Ultima™ Material	
Visco Knee Skin with Straps	4171*
- Closed Patella - Ultima™ Material	
Visco Knee Skin with Straps	6171*
- Closed Patella - Stratus™ Material	

Note: To indicate size, replace \* with corresponding number from the sizing chart on the back cover.

# Q Baby<sup>™</sup>

Compressive strap reduces force on the patellar tendon.

#### **FEATURES**

- Silicone pad may be placed over patellar tendon or tibial tuberosity
- · Fits left or right

#### **INDICATIONS**

- Patellar Tendonitis
- Osgood-Schlatter Disease
- Patellofemoral Knee Pain

Product	Product #
Q Baby Patellar Tendon Strap	5165*

Note: To indicate size, replace \* with corresponding number from the sizing chart on the back cover.

# Premium J™

Closed patella design increases surface contact in the patellafemoral joint.

#### **FEATURES**

- Adjustable, lateral patellar tension
- Specify left or right
- Available with a hinge (see pg. 10)

#### **INDICATIONS**

- Anterior Knee Pain
- · Lateral Tracking Patella
- Patellar Subluxation
- Patellofemoral Knee Pain

Product	Product #
Premium J - Stratus™ Material - Left	6163*L
Premium J - Stratus <sup>™</sup> Material - Right	6163*R



Pull-on Only

## Patella Stabilizer

Stabilizes the patella and provides general support.

#### **FEATURES**

- Foam buttress ring around patella for stability
- · Open patella
- · Fits left or right
- Available with a hinge for increased stability (see pg. 10)

#### **INDICATIONS**

- Patella Stabilization
- Knee Sprains and Strains
- Knee Joint Swelling
- Meniscus Tear
- Osteoarthritis
- Patellar Tendonitis

Product	Product #
Patella Stabilizer - Open Patella - Ultima	a™ Material 4125*
Patella Stabilizer - Open Patella - Stratus	s <sup>™</sup> Material 6125*

Note: To indicate size, replace \* with corresponding number from the sizing chart on the back cover.

# Standard Knee Skin<sup>™</sup>

Compression for minor knee pain and instability.

#### **FEATURES**

- Available in open or closed patella
- Fits left or right
- Available with straps

#### **INDICATIONS**

- Mild Knee Sprains and Strains
- Knee Joint Swelling
- Meniscus Tear
- Osteoarthritis

Product	Product #
Standard Knee Skin - Open Patella	4110*O
- Ultima™ Material	
Standard Knee Skin - Closed Patella	4110*C
- Ultima™ Material	
Standard Knee Skin with Straps	4111*0
- Open Patella - Ultima™ Material	
Standard Knee Skin with Straps	4111*C
- Closed Patella - Ultima™ Material	

Note: To indicate size, replace \* with corresponding number from the sizing chart on the back cover.



## Under Skin<sup>™</sup>

Provides compression and prevents skin irritation while wearing a rigid brace.

#### **FEATURES**

- Thin, breathable compressive material
- 20" in length
- · Limits migration of a rigid brace
- · Latex-free and hypo-allergenic

#### **INDICATIONS**

- Protects Skin from Rigid Braces
- Post-Workout Recovery
- Muscle Injury
- · Mild Knee Swelling

Product	Product #
Underskin	4140*

Note: To indicate size, replace \* with corresponding number from the sizing chart on the back cover.



No Straps

With Straps

## Calf Skin<sup>™</sup>

#### **FEATURES**

- Controls muscle movement
- Reduces muscle vibration

#### **INDICATIONS**

- Contusion
- Shinsplints
- · Calf/Shin Pain

Product	Product #
Calf Skin - I. Iltima™ Material	4630*



## Vector™

Medium-level stabilization for muscle strains and back pain.

#### FEATURES

- Targeted lumbar support
- · Hot/Cold gel pad
- Easy-to-pull, bilateral cinching system
- · Lower profile anterior panel



#### **INDICATIONS**

- Herniated Discs
- Degenerative Disc Disease
- Lumbar Sprains and Strains

Product	Product#
Vector Lumbar Support	3720*

Note: To indicate size, replace \* with corresponding number from the sizing chart on the back cover.



## **Optional** Accessories



Flexible Support



# Back Skin<sup>™</sup>

Compression for muscle strains and minor back pain.

#### **FEATURES**

- Narrow cut in front
- Removable double-cinch belt
- Optional lumbar accessories

#### **INDICATIONS**

- Lumbar Strains
- · Lumbosacral Strains
- Lumbago
- · Lower Back Pain

Product	Product #
Back Skin	4710*
Back Skin with Flexible Support	4734*
Back Skin with Lumbar Foam Pad	4720*
Back Skin with Flexible Support &	4724*
Lumbar Pad	
SI Belt	4758*

Note: To indicate size, replace \* with corresponding number from the sizing chart on the back cover.



## Compression Shorts

#### **FEATURES**

- · Highly Compressive
- Prevents muscle oscillation
- Increases proprioception

#### **INDICATIONS**

- Quad Contusion
- Hip Flexor Strain
- Ischiocapsular/Iliofemoral/ Thigh Strain



#### Product# Compression Shorts - Ultima™ Material 4620\* 4621\* Compression Shorts with Groin Wrap -Ultima™ Material

Note: To indicate size, replace \* with corresponding number from the sizing chart on the back cover.



# Thigh Skin™

#### **FEATURES**

- Gives optimal compression to quads and hamstrings
- Optional Cinch Strap for guad and hamstring injuries

#### **INDICATIONS**

- Quad Contusion
- · Thigh Strain



Product	Product #
Thigh Skin - Ultima™ Material	4610*
Thigh Skin with Cinch Strap - Ultima™ M	aterial 4611*



## Boomerang

#### **FEATURES**

- Universal size
- Easy attachment
- · Soft fleece material against skin
- Fits left or right

#### **INDICATIONS**

- General Wrist Support
- Carpal Tunnel Syndrome
- · Arthritis · Wrist Sprain
- · Wrist Pain · Bursitis

Product	Product #
Boomerang Wrist Wrap	54200



## Standard Elbow Skin<sup>™</sup>

#### **FEATURES**

- Full elbow coverage
- Comfortable, low-profile fit

#### **INDICATIONS**

- Bursitis
- Joint Swelling
- Biceps/Triceps Strains

Product	Product #
Standard Elbow Skin - Ultima™ Material	4210*

Note: To indicate size, replace \* with corresponding number from the sizing chart on the back cover.



# Tennis Elbow Skin™

#### **FEATURES**

- · Silicone pad fits over painful area
- Full elbow coverage

#### **INDICATIONS**

- Medial Epicondylitis (Golfer's Elbow)
- Lateral Epicondylitis (Tennis Elbow)

Product	Product #
Tennis Elbow Skin - Ultima™ Material	4220*
Natar Talia di anta dia manda a * ith	

Note: To indicate size, replace \* with corresponding number from the sizing chart on the back cover.



# Tennis Elbow Band™

#### **FEATURES**

- Silicone pad fits over painful area
- Fits comfortably below elbow

#### **INDICATIONS**

- Medial Epicondylitis (Golfer's Elbow)
- Lateral Epicondylitis (Tennis Elbow)

Product	Product #
Tennis Elbow Band - Ultima™ Material	4250*
Tennis Elbow Band - 2S™ Material	5250*

Note: To indicate size, replace \* with corresponding number from the sizing chart on the back cover.



## Standard Ankle Skin

#### **FEATURES**

- Compressive base sleeve improves proprioception
- · Easily fits into most shoes
- Allows normal range of motion

Product	Product #
Standard Ankle Skin	5350*
Standard Ankle Skin with figure 8 Wrap	5351*

Note: To indicate size, replace \* with corresponding number from the sizing chart on the back cover.




## Visco Ankle

#### **FEATURES**

- Compressive base sleeve improves proprioception
- Includes medial and lateral silicone gel pads
- Easily fits into most shoes
- Allows normal range of motion

Product	Product #
Visco Ankle Skin - Left	5340*L
Visco Ankle Skin - Right	5340*R
V.A.S. w/Figure 8 Compression Wrap - Lef	t 5341*L
VAS w/Figure 8 Compression Wrap - Rig	ht 5341*R





**COMFORTABLE COMPRESSION** 

Powerful Lycra provides a secure fit and perfect compression.



**REFLECTIVE LOGOS & PIPING** 

Stay safe and visible with these strategically placed reflective details.



**EASY STORAGE** 

Large zipper pocket in the back with room for your phone and snacks.

Check out the new colors at: BioSkinFit.com

Page # 6



Q Lok™ ROM
Closed Patella
6166\* - Stratus Material,
Pull-on
5166\* - Front Closure
HCPCS: L1833



CrossFire<sup>™</sup> ROM Closed Patella 6180\* - Stratus Material, Pull-on 5180\* - Front Closure HCPCS: L1833 Page # 6



Gladiator™

Open Patella

4134\* - Ultima Material,
Pull-on

5136\* - Front Closure

6136\* - Stratus Material,
Pull-on

HCPCS: L1833

Page # 5



Hinged Knee Skin™ Open Patella 4130\* - Ultima Material, Pull-on 5133\* - Front Closure 6130\* - Stratus Material, Pull-on

HCPCS: L1820 Page # 7



Closed Patella 6164\* - Stratus Material, Pull-on 5164\* - Front Closure HCPCS: L1820 Page # 9



Visco Knee Skin™ w/ Conforma™ Hinge Closed Patella 6137\* - Stratus Material, Pull-on HCPCS: L1820 Page # 10



Q Brace™
w/Conforma™
Hinge
Open Patella
4136\* - Ultima Material,
Pull-on
5138\* - Front Closure
HCPCS: L1820
Page # 10



Patella Stabilizer™
W/Conforma™
Hinge
Open Patella
4132\* - Ultima Material,
Pull-on
6132\* - Stratus Material,
Pull-on
HCPCS: L1820
Page # 10



w/Conforma™
Hinge
Closed Patella
6138\* - Stratus Material,
Pull-on
HCPCS: L1820
Page # 10

Premium J™



Q Brace<sup>™</sup>
Open Patella
4161\* - Ultima Material
5163\* - Front Closure
HCPCS: A4466
Page # 22



Visco Knee Skin™
Closed Patella
4170\* - Ultima Material,
No Straps
4171\* - Ultima Material,
With Straps
HCPCS: A4466
Page # 22



Q Baby<sup>™</sup> 5165\* HCPCS: A4466 Page # 22



Premium J™
Closed Patella
6163\* - Stratus Material,
Pull-on
HCPCS: A4466
Page # 22



Patella Stabilizer Open Patella 4125\* - Ultima Material 6125\* - Stratus Material HCPCS: A4466 Page # 23



Standard
Knee Skin™
Open or Closed
Patella
4110\* - Ultima Material,
No Straps
4111\* - Ultima Material,
With Straps

HCPCS: A4466 Page # 23



Ten-7™ 5146\* Page # 21



Under Skin™ Closed Patella 4140\* - Ultima Material HCPCS: A4466 Page # 23



Baja<sup>™</sup> 627 3750\* Includes 2 Hot/ Cold Gel Pads HCPCS: L0642 Page # 11



Baja<sup>™</sup> 631 3770\* Includes 2 Hot/ Cold Gel Pads HCPCS: L0648 Page # 11



Baja<sup>™</sup> 637 3760\* Includes 2 Hot/ Cold Gel Pads HCPCS: L0650 Page # 11



Back Flex<sup>™</sup> 3730\* Includes 2 Hot/ Cold Gel Pads HCPCS: L0642 Page # 12



Back Plus™ 3740\* Includes 2 Hot/ Cold Gel Pads HCPCS: L0648 Page # 12



Back Skin™ 4710\* Standard 4720\* w/ Foam

**Lumbar Pad** 4734\* w/ Flexible Lumbar Support

4724\* w/ Foam Pad & Flexible Support

Page # 24



Vector™ 3720\* Includes 1 Hot/ Cold Gel Pad Page # 24



TriLok<sup>™</sup> 5360\* HCPCS: L1902 Page # 15



Visco Ankle Skin<sup>1</sup>

5340\*

5341\* W/Compression Wrap

HCPCS: A4466 Page # 25



Standard Ankle Skin<sup>®</sup>

5350\*

5351\* W/Compression Wrap

HCPCS: A4466 Page # 25



AFTR DC" 5386\* Available with **Compression Wrap** HCPCS: L1902 Page # 17



AFTR ™ 5380\* Available with **Compression Wrap** HCPCS: L1902 Page # 17



AFTR w/ Gel™ 5383\* Available with **Compression Wrap** HCPCS: L1902 Page # 17



Hallux Control Strap<sup>™</sup> 9534\*T Ten Pack CPT 29550 Page # 18



Weil Osteotomy Strap<sup>®</sup> 95372T Ten Pack CPT 29550



Calibrated Midfoot Compression Wrap<sup>'</sup> 95362T Ten Pack CPT 29540 Page # 18



DP2<sup>™</sup> 6 1/2" Cock-Up Wrist 5423\* HCPCS: L3908 Page # 13



DP3<sup>™</sup> 8 1/2" Cock-Up Wrist 5424\* HCPCS: L3908 Page # 13

Page # 18



Wrist/Thumb Spica 5452\* HCPCS: L3807 Page # 13



Thumb Spica 5450\* HCPCS: L3923 Page # 13



Boomerang<sup>"</sup> 54200 HCPCS: A4466 Page # 25



Standard Elbow Skin™ 4210\* - Ultima Material HCPCS: A4466 Page # 25



Tennis Elbow Skin<sup>®</sup> 4220\* - Ultima Material HCPCS: A4466 Page # 25



Tennis Elbow Band™ 4250\* - Ultima Material 5250\* - 2S Material HCPCS: A4466 Page # 25



Calf Skin™ 4630\* - Ultima Material HCPCS: A4466 Page # 23



Thigh Skin™ 4610\* - Ultima Material 4611\* - Ultima Material, with Cinch Strap HCPCS: A4466 Page # 24



Compression Shorts 4620\* - Ultima Material 4621\* - Ultima Material, with Groin Wrap

HCPCS: A4466 Page # 24




### **Application Instructions**




See product page on-line.

#### Knee



	Knee Products				
		Circumferen	ce with knee fu	lly extended	
Size	*PN	Incl	hes	Centimeters	
		mid – patella	7" above	mid – patella	18cm above
XS	0	12 – 13	16 – 18	31 – 33	41 – 46
S	1	13 – 14	18 – 20	33 – 36	46 – 51
М	2	14 – 16	20 – 22	36 – 41	51 – 56
L	3	16 – 17	22 – 25	41 – 43	56 – 64
XL	4	17 – 19	25 – 28	43 – 46	64 – 71
2XL	5	19 – 21	28 – 32	46 – 54	71 – 81
3XL	6	21 – 23	32 – 36	53 – 58	81 – 91
4XL	7	23 – 25	36 – 40	58 – 64	91 – 102
5XL	8	25 – 27	40 – 44	64 – 69	102 – 112




	Ten-7				
	Circumference of upper and lower thigh				
Size	Size *PN Inches Centimeters			neters	
		Upper thigh	Lower thigh	Upper thigh	Lower thigh
1	1	18.5 – 19.5	14 – 15.5	47 – 50	36 – 39
2	2	19.5 – 21.5	14.5 – 16.5	50 – 55	37 – 42
3	3	21.5 – 23.5	15.5 – 18	55 – 60	39 – 46
4	4	23.5 – 26	17 – 19.5	60 – 66	43 – 50
5	5	26 – 28.5	18.5 – 21	66 – 72	47 – 53
6	6	28.5 – 31	20 – 22.5	72 – 79	51 – 57



	Q Baby			
	Circumference under patella with knee fully extended			
Size *PN Inches Centimeters				
XS - M	2	12 – 16	31 – 41	
L-XXL	5	16 – 21	41 – 54	

### Thigh




Mid — thigh circumference				
Size *PN Inches		Inches	Centimeters	
XS	0	16 – 18	41 – 46	
S	1	18 – 20	46 – 51	
М	2	20 – 22	51 – 56	
L	3	22 – 25	56 – 64	
XL	4	25 – 28	64 – 71	
2XL	5	28 – 32	71 – 81	

#### Calf



Largest part of calf circumference			
Size *PN Inches			Centimeters
S	1	11 – 12	28 – 30
M	2	12 – 14	30 – 35
L	3	14 – 16	35 – 41
XL	4	16 – 18	41 – 46
2XL	5	18 – 20	46 – 51

### Spine & Shorts



Circumference between hip and top of pelvis.					
Size	*PN	Inches	Centimeters		
Baja™ & Vector™					
S	1	26 – 30	66 – 76		
М	2	30 – 36	76 – 91		
L	3	36 – 44	91 – 112		
XL	4	44 – 54	112 – 137		
UNIV	4U	26 – 54	66 – 137		
Back Flex & Back Plus					
1	1	28 – 34	71 – 86		
3	3	34 – 44	86 – 112		
5	5	44 – 56	112 – 142		
Back Skin					
XS	0	24 – 26	61 – 66		
S	1	26 – 30	66 – 76		
М	2	30 – 34	76 – 86		
L	3	34 – 38	86 – 97		
XL	4	38 – 42	97 – 107		
2XL	5	42 – 46	107 – 117		
Compression Shorts					
S	1	28 – 32	71 – 81		
М	2	32 – 36	81 – 92		
L	3	36 – 40	92 – 102		
XL	4	40 – 44	102 – 112		
2XL	5	44 – 48	112 – 122		

#### **Elbow**



Circumference of largest part of forearm					
Size	*PN	Inches	Centimeters		
XS	0	8 – 9	20 – 23		
S	1	9 – 10	23 – 26		
М	2	10 – 11	26 – 28		
L	3	11 – 12	28 – 31		
XL	4	12 – 13	31 – 33		
2XL	5	13 – 14	33 – 36		

Call for custom products and extended sizing.

### Indicate Sizing When Ordering



To indicate size, replace \* with corresponding \*PN (Part Number) from the sizing chart.


**Continued on Back Cover** 






## SIZING CHART

#### Foot & Ankle





#### Wrist

000	Circumference of extended wrist			
	Size	*PN	Inches	Centimeters
	XS - S	1	4 – 6	10 – 15
- M-	M-L	3	6-8	15 – 20
//// \	XL-2XL	5	8 – 10	20 – 25

٥٩٥	Thumb Spica			
	Length of thumb			
\	Size	*PN	Inches	Centimeters
	XS - M	2	3.25 – 4.25	8 – 11
·	L-2XL	5	4.25 – 5.25	11 – 14

Call for custom products and extended sizing.

# ON-LINE ORDERING

USA: <u>Pro.BioSkin.com</u>
UK: <u>BioSkinUK.co.uk</u>

Easy, Fast, Secure, 24/7

- Place Orders
- Pay Invoices
- Application Videos
- Product Reviews



USA: Orders@BioSkin.com
UK: UK@BioSkin.com

Email orders 24/7

Use Subject Line: "Bio Skin Order"



USA: 1-800-541-2455

Monday - Friday, 7 am - 5 pm Pacific time.

UK: +44 (0) 845 528 0276

Monday - Friday, 08.30 - 17.30



USA: 1-541-482-2119

UK: +44 (0) 148 342 1333

Indicate Sizing When Ordering



To indicate size, replace \* with corresponding \*PN (Part Number) from the sizing chart.

**More Sizing on Previous Page** 

