


TherAlign Dragonfly TLSO • 2018 Sizing Chart


STEP 1: Record measurements

Garment Subsystem



Dragonfly Brace

Spine Length C-7 to Coccyx



STEP 2: Size & Select Garment Subsystem

SIZE GROUP	WEIGHT (LBS)	STA (IN.)	ATK (IN.)	GARMENT SUBSYSTEM
Infant	15 – 25	8 – 10	7 - 10	☐ 110
Pre-School	26 – 45	11 - 13	11 - 12	☐ 210
Pediatric	46 – 65	14 - 15	13 - 18	☐ 310
Junior	66 – 120	16 - 17	19 - 21	☐ 410
Female Med.	121 – 170	18 – 20	22 - 24	☐ 610
Female Large	171 – 200	21 - 22	25 - 27	☐ 620
Male Med.	121 – 170	18 – 19	22 - 24	☐ 710
Male Large	171 – 190	20 - 22	25 - 27	☐ 720

STEP 3: Size & Select Dragonfly Brace

SPINE (IN)	MODEL
10 - 12	☐ DFLY30
13 - 14	☐ DFLY33
15 - 16	☐ DFLY40
17 – 19	☐ DFLY48
20 - 21	☐ DFLY53

STEP 4: Size/Select Rib and Pelvis Strapping

ADULT: 2-INCH STRAPPING

CHEST GIRTH (IN)	RIB STRAP	HIP GIRTH (IN)	PELVIS STRAP
23 - 30	☐ PS-204	17 - 22	☐ PS-204
31 - 40	☐ PS-206	23 - 30	☐ PS-206
41 - 54	☐ PS-210	31 - 41	☐ PS-210
55 - 62	☐ PS-214	42 - 48	☐ PS-214

PEDIATRIC: 1-INCH STRAPPING

CHEST GIRTH (IN)	RIB STRAP	HIP GIRTH (IN)	PELVIS STRAP
16-24	☐ PS-103	15-19	☐ PS-103
25-28	☐ PS-105	20-25	☐ PS-105
29-33	☐ PS-107	26-30	☐ PS-107

Rib & Pelvis Strapping

