

DUAL ACTION Articulating Ankle Foot Orthosis



**Posterior ankle strap is shown inside upright for photo. When fitted properly, the strap is outside posterior upright bar.

Features & Uses

- ▶ The Dual Action Orthosis offers patients the same outstanding features and uses of our APU® (Adjustable Posterior Upright) Orthosis and the EV™ (Equino-Varus) Orthosis in one, with + or - 20° adjustments in Dorsiflexion and + or - 40° in Plantarflexion and Varus/Valgus adjustability as well.
- ▶ Offers measurable adjustability features for the ambulating patient
- ▶ Proper support and heel suspension for the recumbent patient.



Suggested L-4396

MADE IN THE
USA

ANATOMICAL
CONCEPTS, INC.