


To determine proper liner sizing, one measurement is required:

**Step #1:** Take circumference measurement 2" (5-6cm) proximal to the distal end of the residual limb.

**Note:** Unless specified, all fabric covered Seal-Pro TF liners come with proximal 7" exposed.  
If custom fabric is desired, please state your description at the time of order.

STEP  
1



STEP 1: Limb Circ.	Choose Part Number
24 - 28 cm	SP-TF-22 - <b>A</b>
28 - 32 cm	SP-TF-26 - <b>A</b>
32 - 36 cm	SP-TF-30 - <b>A</b>
36 - 40 cm	SP-TF-34 - <b>A</b>
40 -44 cm	SP-TF-38 - <b>A</b>
44 - 48 cm	SP-TF-42 - <b>A</b>
48 - 52 cm	SP-TF-46 - <b>A</b>

STEP 2: Custom Textile Requests	Insert Code for "A"
Cushion NO Textile	* Leave Blank *
Cushion w/ Textile Prox. 7" Exposed	(F)
Cushion w/ Custom Textile to Spec.	(F) Custom
Locking w/ Textile Prox. 7" Exposed	(L)
Locking w/ Custom Textile to Spec.	(L) Custom

**EXAMPLE:** SP-TF-34 - **L Custom**

Limb Circ. Of 37" and a Lock with Custom Textile.