



APRL Child Flexion Wrist


HOSMER

- Installs between terminal device and wrist allowing manual positioning of terminal device
- Flexes at 0 and 36 degrees
- 1-1/16" (2.7cm) long, 1.5oz (42g)
- Recommended for use with WE-100N
- Standard 1/2-20 thread and stub

52176 APRL Child Wrist Flexion Unit WF-50

Includes:

- A. 52177 Body
- B. 52179 Rocker
- C. 52182 Friction Plug
- D. 52185 Set Screw
- E. 52183 Screw
- F. 52181 Plunger
- G. 52180 Spring
- H. 72056 Set Screw



Sierra Wrist Flex Unit

HOSMER

- Lever selects three positions, 0, 25 or 50 degrees of flexion
- 1-7/16" (3.6cm) long, 4oz (113g)
- Standard 1/2-20 thread and stud

52371 Sierra Wrist Flex Unit 18-00

Includes:

- A. 52380 Base
- B. 52375 Lever
- C. 52466 Pin
- D. 54954 Rocker
- E. 52372 Rocker Bushing
- F. 52379 Cover
- G. 52301 Screw
- H. 52790 Pin
- I. 52378 Plug
- J. 52300 Set Screw
- K. 52376 Compression Spring

