



2017 Catalog (Update)

**FOOTWEAR FOR DIABETICS, ORTHOPEDICS,
PEDORTHICS & PROSTHETICS
CUSTOM MADE SHOES, BOOTS & INSOLES
FOR HARD-TO-FIT NEEDS**



Apis Footwear Co. Tel 1.888.937.2747 Fax 1.888.990.2245
2239 Tyler Avenue, South El Monte, CA 91733
13026 Forest Centre Court, Louisville, KY 40223
www.apisfootwear.com. email: sales@apisfootwer.com

A PROUD MEMBER OF AOPA, AAOP & PFA

“Wide or narrow, big or small,
Apis has a shoe for everyone”.

Toddlers: 4, 4.5 - 11 in M, W, XW

Youth: 11, 11.5 - 6 in M, W, XW

Ladies: 4.5 - 15 in AA, B, 3E, 5E, 7E, 10E, 14E

Men: 5 - 20 in D, 2E, 4E, 6E, 9E, 10E, 14E

Whether for walking, casual dress, athletic
or medical therapeutic footwear

We have what you need,
off-the-shelf or custom made.
Shoes or boots, foot orthotics,
functional or accommodative.

Women's Premiere Collection

Medicare approved, feature seamless linings, true width and added depths, better accommodation for orthotics, AFO & internal modifications, better support & stability.

Features:

Genuine leather upper with seamless lining

Added width at toe & extra support at the ankle

Added depth - 8mm or 5/16" deeper than conventional straight last shoe

Better accommodation for orthotics, AFO & internal modification

Long medial counter to provide support & stability

Removable inserts for flexible fitting

Widths: B/D(M/W), 3E/(WW), 5E/(XW) **Sizes:** 5, 5.5 - 10.5, 11



9205-Ruby



9205-Brown



9205-L-Black



9205-Black

(Elastic Vamp - ideal for Bunion, Hammer Toes, mild to moderate swelling situation)



9200-Black



308 (A5512)

Women's Premiere Collection



9212 Taupe



9212 Black



9214 Black Lycra
(Elastic Vamp - accommodate
cross-over toes or Hammer toes)

Sizes & Widths:

B/D(M/W): 5 - 10.5, 11

3E(WW): 5 - 10.5, 11

5E(XW): 5 - 10.5, 11

Genuine leather uppers
True widths with wide toe box
Extended medial & lateral counter
Reinforced support & stability
Seamless fabric linings



9215 Black 5" Boots

Genuine leather boots with elastic
laces and single strap,
Diagonal zipper for easy entry
Extended counter for reinforced stability

Breathable Mesh Walking Shoes -Widths: D, 3E, 5E Sizes: 5, 5.5 - 11.5, 12

(Breathable Mesh Vamp - accommodate
oot orthotics, AFOs, prosthetics and internal
modifications)



3310 Black Mesh Walking

Men's Premiere Collection

Genuine leather upper
 High density injected EVA outsole
 Symmetrical cut-line for easy sole modifications
 Removable inserts for flexible fitting
 Added depth - 8mm or 5/16" deeper than conventional last
 Added width at the toe & extra support at the ankle
 Better accommodation for orthotics, AFO & modifications
 Seamless linings to decrease friction
 Feather-light - weighs only 8 ounces per shoe



Men's Casual Dress

708-L Casual Dress Lace



888-V Casual Dress Velcro



9602-Brown Casual Dress Velcro



9602-Black Casual Dress Velcro Black Lycra

Sizes & Widths:
 D/2E(M/W): 7 - 12.5, 13 - 15
 4E(XW): 7 - 12.5, 13 - 15
 6E(XXW): 7 - 12.5, 13 - 15

High Instep/ Swelling	Pronated/ Flat Foot	Bunions	Hammer Toes	Partial Foot	AFO	Weight	Closure
Moderate	Moderate	Moderate	Moderate	Yes	Yes	8 oz.	Velcro/ Lace

Men's Premiere Collection

Men's Athletic Walking

9701-1L/Black



9701-1V/Black



9701-3L/White/Silver



9701-3V/White/Silver



9701-5L/Grey/Silver



9701-5V/Grey/Silver



Widths: D/2E(M/W), 4E(XW), 6E(XXW)

Sizes: 7, 7.5 - 12.5, 13 - 15

Standard American Last (FIA) , match Ritzstick or Brannock sizing

Genuine leather upper and fabric linings

High Density solid EVA outsole

Seamless linings, removable insoles

Added depth - 8mm or 5/16" deeper than conventional shoes

Wider toe boxes and reinforced counter

Accommodation of foot orthotics, AFOs & internal modifications

PDAC A5500 approval

Men's Premiere Collection

Men's Casual Boots



9605 Black



9605 Wax Tan



9606 Lace Black



9606-V Black Hook & Loop
w/zippers

Widths: D/2E(M/W), 4E(XW), 6E(XXW)
Sizes: 7, 7.5 - 12.5, 13 - 15

Genuine leather upper
High density injected EVA outsole
Symmetrical cut-line for easy sole modifications
Removable inserts for flexible fitting
Added depth - 8mm or 5/16" deeper than conventional last
Added width at the toe & extra support at the ankle
Better accommodation for orthotics, AFO & internal modifications
Seamless linings to decrease friction
Feather-light - weighs only 8 ounces per shoe

High Instep/ Swelling	Pronated/ Flat Foot	Bunions	Hammer Toes	Partial Foot	AFO	Weight	Closure
Moderate	Moderate	Moder- ate	Moderate	Yes	Yes	8 oz.	Velcro/ Lace

Women's Comfort Collection

Casual Dress



9202

Sizes & Widths:
 AA: 4.5, 5 - 11, 12
 B: 4.5 - 10.5, 11 - 15
 D: 4.5 - 10.5, 11 - 15
 3E: 4.5 - 10.5, 11 - 15
 5E: 5 - 10.5, 11 - 15
 7E: 5 - 10.5, 11 - 15

Genuine leather
 Double depth combination last
 3 layers of removable inserts
 Extended firm heel counter
 Solid EVA outsole
 Great for AFOs and braces



9301 Black



9302 Black



9106 Black

Oblique Toe
 Supra depth/Triple Depth



9108 Black

Breathable Mesh Walking Shoes -Widths: D, 3E, 5E Sizes: 5, 5.5 - 11.5, 12

9306 - Black



9315 - Black



High Instep/ Swelling	Pronated/ Flat Foot	Bunions	Hammer Toes	Partial Foot	AFO	Weight	Closure
Severe	Severe	Severe	Severe	Yes	Yes	10 oz.	Velcro/ Lace

Men's Comfort Collection

Men's Casual 801/802 Supra-Depth

801 Black



801 Brown



802 Black



Sizes & Widths:
D: 5 - 11.5, 12 - 17
4E: 6 - 11.5, 12 - 17
6E: 7 - 11.5, 12 - 17
9E: 7 - 11.5, 12 - 17

802 Brown



7021 Black



Triple / Supra depth
Reinforced Counter and toe box
3 layers of removable inserts
High density solid EVA outsole
Genuine soft leather uppers & lining

Sizes & Widths
D: 5 - 11.5, 12 - 18
2E: 5 - 11.5, 12 - 18
4E: 6 - 11.5, 12 - 20
6E: 7 - 11.5, 12 - 20

Breathable Mesh Walking Shoes - Widths: D, 4E, 6E, & 9E
Sizes : 7, 7.5 - 11.5, 12 - 17

9704 Black



9713 Boots Black



Men's Comfort Collection

Men's Casual Dress 9501 / 9502 / 502 Double-Depth

Sizes & Widths:

D: 5 - 11.5, 12 - 18

2E: 5 - 11.5, 12 - 18

4E: 6 - 11.5, 12 - 20

6E: 6 - 11.5, 12 - 20

9E: 6 - 11.5, 12 - 20

Genuine leather upper

High density injected solid EVA outsole

Removable inserts for flexible fitting

Added width at the toe & extra support at the ankle

Accomodate orthotics, AFO & internal modifications

Seamless linings to decrease friction

Feather-light - weighs only 8 ounces per shoe



9501 Black



9502 Black Slap-Over strap



502 Black



502 Brown



9921 Black



9921 Tan

Genuine leather upper & lining
High density solid EVA outsole
Oil-resistant outsoles
Removable inserts for flexible fitting
Double depth & extended width

Sizes & Widths:

D: 6 - 11.5, 12 - 17

4E: 6 - 11.5, 12 - 20

6E: 7 - 11.5, 12 - 20

Men's Comfort Collection

9703-2L Hiking Boots

Suede and mesh Vamp
Solid Injected TPR outsole for durability
Removable insoles and reinforced counter
Added depth and wider toe box
Better accommodation of orthotics,AFOs
and modifications.
Ideal for ankle support and protection

Widths & Sizes

D: 7 - 11.5, 12 - 17

4E: 7 - 11.5, 12 - 17

6E: 7 - 11.5, 12 - 17

9E: 7 - 11.5, 12 - 17



AM5605 Lace

6" Adventure Boots

Widths: D/W,4E/XW, 6E/XXW

Sizes: 7, 7.5 - 12.5, 13 -15



AM5615 Straps

5"Chukka Boots - Duo -Depth

Widths: D,4E, 6E

Sizes: 5, 5.5 - 11.5, 12 - 20

Quality Genuine leather
Extra-Depth Boots
Three layers of removable insoles
Solid EVA outsole
Reinforced counter



9951 Black



9951 Tan

Men's Comfort Collection

- Widest Athletic Shoes With 9E width

9702-1V



9702-3V



9702-1L



9702-3L



Widths & Sizes

D: 7 - 11.5, 12 - 17

4E: 7 - 11.5, 12 - 17

6E: 7 - 11.5, 12 - 17

9E: 7 - 11.5, 12 - 17

Genuine leather uppers and breathable mesh vamp

Solid TPR injected outsole, reinforced counter,

Removable insoles, double depth shoe with wider toe box

Ideal for foot orthotics, AFO's, prosthetic foot, braces and modifications

9708-1L



9708-2L



9708-3L



Widths & Sizes

D: 7 - 11.5, 12 - 17

4E: 7 - 11.5, 12 - 17

6E: 7 - 11.5, 12 - 17

9E: 7 - 11.5, 12 - 17

Suede leather with breathable mesh vamp

width extends to 9E, wider toe box

Wide opening for easy entry of
foot orthotic devices

Men's Comfort Collection



504 - Black



503 - Black

504 8"- Boots with zipper & laces

Genuine leather upper
New improved medial zipper on 504
Supra depth with wide toe box
Three layers of removable insoles
Reinforced counter for stability

503 8"- Boots with double straps

Widths & Sizes

D: 5 - 11.5, 12 - 17
2E: 5 - 11.5, 12 - 17
4E: 6 - 11.5, 12 - 17
6E: 7 - 11.5, 12 - 17
9E: 7 - 11.5, 12 - 17



6506 - Utility Protection 8" Boots (PDAC & OSHA)

Genuine leather, reinforced counter
OSHA approved Composite Toe cap.
Double depth with wide toe box.
Widths: D, 4E, 6E Sizes: 6, 6.5 - 11.5, 12 - 15

Men's Oil-resistant / Anti-slip Utility Shoes

--- Ideal for garage, restaurant, hotel and sanitary industries

2011

2012

2013



Features:


Utility shoe collection
Genuine leather upper
oil / slip resistant rubber soles
removable insoles
Added depths
Three true widths

Sizes & Widths:

D: 7, 7.5-11.5, 12-13
4E: 7, 7.5-11.5, 12-13
6E: 7, 7.5-11.5, 12-13



FOR HARD-TO-FIT PATIENTS - MEDICAL COLLECTION



Edema	9301-E, 9301-X, 502-E, 502-X, 628-E, 728-E
Bunion/Hammer Toe/Crossover	608, 618, 628-E, 708, 718, 728-E
Charcot/CMT	9301-C, 502-C, 503-C, 504-C
Surgical opening	511, 9226

Women's Medical Collection

Medical Collection of shoes are designed specifically for patients who need larger, wider shoes. These shoes offer variety of closures and material combinations to accommodate hard-to-fit foot, such as, charcot, edema, hammer toes or foot trauma.

628-E Bunion / Edema for woman



Cut-out flex area to release bunion pressure
Padded collar reduces friction on the heel
Elastic vamp accommodates hammer toe or swelling

Sizes: 4.5, 5 - 10.5, 11 - 15
Widths: B, D, 3E, 5E, 10E, 14E

608, 618 Bunion / Bunionette

Widths: B, D, 3E, 5E

Sizes: 4.5 - 10.5, 11 - 15



608 Black



618 Black

Cut-out flex area to release bunion pressure
Padded collar reduces friction on the heel
Elastic vamp accommodates hammer toe or swelling

9226 - Surgical Opening



Widths: D, 3E, 5E
Sizes: 5, 5.5 - 11.5, 12 - 15

Genuine leather
Surgical opening for easy entry
Built-in solid EVA rocker sole
Perfect for braces, orthotics and modifications

For Charcot, Edema, Hammer Toes & Pes Planus

9301-C Charcot Black



Widths: B, D, 3E, 5E, 7E
Sizes: 4.5 - 11.5, 12 - 15

Genuine leather upper
Seamless linings
Solid EVA outsole
Mis-mates 9301 or 9301 - X
Extended width and volume at arch



The Accommodator 9301- X / 9301- E

9301-E features expandable elastic foam cushion collar for sensitive ankle/instep severe edema, comes in 10E and 14E, the widest available width of its kind. This collar can be added to other styles in the catalog as modification order.



9301-E Black



9301-X Black

Widths: 10E (W), 14E(EW)
Sizes: 5, 5.5 - 11.5, 12 - 15

Genuine leather upper
Added Supra Depth
Wide entry opening
Built-in Volume / Space
Mis-mate 9301 or 9301-C
Fits up to 25" ankle girth

Women's Medical Collection

- Wound Care Shoe Selection

Expandable Lycra vamps with padded collar
Accommodate severe edema, sensitive skin & R.A
Post operations care and healing protection

638 (For Woman)



Widths & Sizes

M: 5 - 11.5, 12 - 15
W: 5 - 11.5, 12 - 15
XW: 5 - 11.5, 12 - 15

639 (For Woman)



Leather straps
Stretchable lycra vamps
Accommodate bunion/bunionette
Hammer toes / overlapping, claw toes, Edema
Solid EVA outsole for easy modification
Three true widths and size up to 15

637 (For Woman)

Stretchable Lycra vamps with adjustable strap



Widths & Sizes

D: 5 - 11.5, 12 - 15
3E: 5 - 11.5, 12 - 15
5E: 5 - 11.5, 12 - 15

Men's Medical Collection

Bunion / Edema Shoes 728-E

Cut out flex area for bunion relief

Padded foam collar reduces friction around the ankle

Flex leather to reduce achilles pressure

Expandable vamps offers tremendous forgiveness to toe pressure & swelling

Multitude of widths and sizes

Designed for severe bunion / edema



Widths & Sizes

D: 5 - 11.5, 12 - 17

4E: 5 - 11.5, 12 - 17

6E: 5 - 11.5, 12 - 17

10E: 5 - 11.5, 12 - 17

14E: 5 - 11.5, 12 - 17

708, 718 Bunion / Bunionette

708 Black



718 Black



Widths & Sizes

D: 5 - 11.5, 12 - 17

4E: 5 - 11.5, 12 - 17

6E: 5 - 11.5, 12 - 17

511 Black - Surgical Opening



Widths & Sizes

D: 5 - 12.5, 13 - 15

3E: 5 - 12.5, 13 - 15

5E: 5 - 12.5, 13 - 15

Genuine leather upper

Double depth

Surgical opening entry

Built - in rocker sole

Perfect for AFOs or foot orthotics

502-C, 503-C, 504-C Charcot last

502-C Black



503-C Black



504-C Black



Size: 5, 5.5 - 11.5, 12-17
Width: D, 2E, 4E, 6E, 9E

Medial Charcot last
Seamless toe box
Solid EVA sole
Mis-mates 502, 502 - X or 502 - E

The Accommodator 502-X / 502-E

Extra supra-depth to accommodate severe edema / swollen feet

502-X Black



502-E Black



Genuine leather upper
Heel counter cushioned with elastic foam
Designed for sensitive ankle / instep severe edema
Offered in 10E & 14E
Fit up to 25" ankle

Size: 5, 5.5 - 11.5, 12-17
Width: 10E, 14E

Men's Medical Collection

For post-op protection, severe Edema, Sensitive skin & R.A



738(For Men)

M: 5 - 11.5, 12 - 15
W: 5 - 11.5, 12 - 15
XW: 5 - 11.5, 12 - 15



739 (For Men)

Stretchable Lycra Vamp
Leather straps
Accommodate bunion / Hammer toe / Edema
Solid EVA outsole for easy modification
Three true widths

Post-op Lycra Vamps with Adjustable strap



737 (For Men)

Widths & Sizes
D: 5 - 11.5, 12 - 15
4E: 5 - 11.5, 12 - 15
6E: 5 - 11.5, 12 - 15

Mt Emey Diabetic Socks

Extra soft bamboo fiber
Non-binding crew socks
No-sagging arch area
No seam design
Loose fit to make sure good circulation
Moisture control to reduce bacteria infection
Available in black and white



Children's Orthopedic Shoes



3301-1L

Features

- Straight-last
- Extra-depth with 4mm
- Removable inserts
- Elongated counter support
- Wide opening for easy entry
- Elastic laces & D-ring hook & loop
- Non-marking soft rubber outsole
- Solid high density EVA midsole

Sizes & Widths

Sizes: Toddler (T): 10-13 (full sizes)

Youth (Y): 1-5, 5.5-6.5

Men(Y): 7-10.5, 11

Widths: M, W, XW

Colors: White, Black, Earth, Blue, Pink



3301-3L



3301-5L



3301-6L



3301-7L

Children's Orthopedic Shoes



3305-1H

Features

- Straight-last
- Extra-depth with 4mm
- Removable inserts
- Elongated counter support
- Wide opening for easy entry
- Elastic laces & D-ring hook & loop
- Non-marking soft rubber outsole
- Solid high density EVA midsole

Sizes & Widths

Sizes: Toddler (T): 10-13 (full sizes)

Youth (Y): 1-5, 5.5-6.5

Men(Y): 7-10.5, 11

Widths: M, W, XW

Colors: White, Black, Earth, Blue, Pink



3305-3H



3305-5H



3305-6H



3305-7H

Children's Orthopedic Boots



2684



(2684
Outsole)

(2681 & 2683
Outsole)

Features

- Straight-last
- Extra-depth with 4mm
- Removable inserts
- Elongated counter support
- Wide opening for easy entry
- Solid high density EVA midsole
- Non-marking soft rubber outsole

Sizes & Widths

Sizes: Toddler (T): 4-13.5 (half sizes)
Youth (Y): 1-6 (half sizes)
Widths: M, W, XW
Colors: White, Yellow, Pink



2681



2683

Children's Semi-Surgical Footwear



2603-V



2603-L

Special Order Category
Non-Stock Items

Features

- Semi-surgical opening
- Adjustable depth
- Removable inserts
- Elongated counter support
- Wider heel & toe box for orthoses or AFO/KAFO
- High abrasion resistant injected rubber sole
- Fray resistant rear foot linings

Sizes & Widths

Sizes: Toddler (T): 4-11, 11.5-13.5
Youth (Y): 1-4, 4.5-6.5
Men (Y): 7-10.5, 11-14
Width: M, W, XW
Color: White



2623
(Navicular strap 5"boots)



2633
(U-type heel entry boots)

Children's Non-slip Footwear



2151

Features

- Oil & slip resistant rubber outsole
- Straight-last
- Extra-depth with 4mm
- Removable inserts
- Elongated counter support
- Wide opening for easy entry

Sizes & Widths

Sizes: Toddler (T): 11-13.5 (half sizes)

Youth (Y): 1-6 (half sizes)

Widths: M, W, XW

Colors: Black, Pink, Blue



2155



2153

(non-stock item)



2152



2154

(non-stock item)



2156

(non-stock item)

Children's Stability Footwear



2571-SL

Features

- Straight-last (SL)
- Extra-depth
- Removable inserts (2571 , 2572. 2574)
- Elongated counter support
- Non-marking, gum rubber outsole
- Toe-In (TI) and Toe-Out (TO) last are also available as non-stock item



2572-SL

Sizes & Widths

Sizes: Toddler (T): 11-13.5 (half sizes)

Youth (Y): 1-6 (half sizes)

Widths: M, W, XW

Colors: Black, Blue, Pink



2574-SL
(non-stock item)



2573-SL
(non-stock item)



1368-SL
(non-stock item)



1369-SL
(non-stock item)



(Ousole)

Toddler Shoes



Widths: Medium / B
Wide (EE)
Sizes: 4 - 11.5



222-3 (White / Pink) Girl

222-3 (White / Navy) Boy

Youth Sneakers



223-3 White/Navy
Boy's



224-3 White/Pink
Girl's



226-3 White/Silver
(unisex - adolescent)



230-1 Black



231-3 White



229-1 Black



229-3 White

Genuine leather upper
Symmetrical straight last
Added depth - 5/16" deeper than conventional last
Reinforced counter for stability and support
removable insoles for flexible fitting
Extra overall volume
Solid EVA outsole for easy modification
Ideal for AFOs, UCBLs and braces

Widths: Medium / B
Wide (EE)
Sizes: 11, 11.5 - 13.5, 1 - 6

Women's Answer2 Orthopedic Shoes

ANSWER²™

Walking Shoes

447-1 Black



Widths: M, W, XW
Sizes: 5, 5.5 - 10.5, 11

447-3 White



448-1 Black



448-3 White



Genuine leather upper
Double depth last
Semi-surgical opening
Reinforced counter
Solid EVA outsole for easy modification
Perfect for UCBLs, AFOs, Foot orthotics & prosthetics

Casual Shoes



445-1 Black



446-1 Black

Widths: M, W, XW
Sizes: 5, 5.5 - 10.5, 11

Athletic Shoes

554-3 White/Navy



554-3V White/Navy



554-1 Black



553-1 Black



558-1 Black



558-3 White



Genuine leather upper
Double depth last
Semi-surgical opening
Reinforced counter
Solid EVA outsole for easy modification
Perfect for UCBLs, AFOs and braces

Sizes: 6 - 12.5, 13 - 15
Width: Medium (D)
Wide (2E)
ExtraWide (4E)

Men's Stability Boots - Broad Entrance

551-1 Black



552-1 Black



Custom Shoe Program



Risk-free program, guaranteed fit

One flat price, any modifications

Any style in the catalogue, over 70 styles available

Heel lift
Fore foot rocker
Bevel heel
SACH heel
High top support
Stabilizer

Rigid sole
Semi-rigid sole
Heavy-duty Rubber outsole
3/4 length steel shank
Full length steel shank

Semi-surgical opening
Surgical opening
Composite toe box
Thermozote lining
Padded tongue or collar
Leather lining
and more...

Custom Foot Orthotics

Accommodative Or Functional

Risk-free program, guaranteed fit



L5000 PDAC Toe Filler
513 (A5513) Custom Molded Diabetic Insoles



Custom Functional Orthotics
523 - Full Length or 3 / 4




UCBL

For Patients With Diabetes And Foot Pain

Choice of base and top cover materials
EVA, thermal-cork, rubber cork, Thermo-zote, PPT, Spenco,
Heat molded over positive foot model of each patient
Complies with new Medicare guidelines

Mt Emey Prefab Orthotics

... Ultimate Comfort & Support



Mt. Emey OTS Foot Orthotics provide the ultimate support and comfort for patients while simultaneously yielding the highest profit margin for practitioners. Our goal is to create a dual satisfying product. It offers built-in arch support and heel pad for heightened support, tri-lam layer cushioning with our own EVA base for unique comfort and durability. Extremely durable and firm shell for protection and corrective support. At a really great price per pair, you can begin to introduce a new line of comfort for both you and your patients.



Inserts/Orthotic Blanks

Pre - fab Accommodative Inserts

308 Insert (men & women)
1/8" thickness Thermozote cover; 1/4" thickness heat-moldable close-cell EVA base comes in white, or blue, or green color.



Reference Table

Accommodative Inserts		
308/92	Lady's	9201,9202,8821,8822,8824,8825
308/91	Lady's	9106,9108(oblique toe)
308/93	Lady's	9301,9302,608,618,628-E, 8935, 8936, 8937, 408,9205,9206,9207,9208,9209,9210,9211, 9212,9213,9203.
308/95	Men's	502,508,509,7021,7022,708,718,888,9501,9502,9504,9505,9901,9912,9921,9923,9925,9928, 9931,9938,9951,503,504,801,802
305	Lady's & Men's	9301-C,9301-R,9301-X,502-C,502-E,502-X, 503-C,504-C

For diabetics only, multiple density insert, direct formed, molded to foot after external heat source of about 230F or higher, total contact with foot, including arch, base layer min. of 1/4" material of Shore A35 durometer or 3/16" material of Shore A 40 (or higher).

Orthotics Blanks for Custom Orthotics



431-E Blanks:

Top: 3/16" Thermozote
Mid-layer: 1/8" PPT/EVA 30 - 35 Shore A
Arch: Semi-oval crescent shape EVA in 35 - 45 Shore A tapered to lateral
Base: 1/4" close cell EVA in 50 - 55 Shore A, tapered to toe 1/8"

432-E Blanks:

Top: 3/16" Thermozote
Mid-layer: 1/8" PPT/EVA 30 - 35 Shore A
Arch: Semi - oval crescent shape EVA in 35 - 45 Shore A tapered to lateral
Base: 1/2" close cell EVA in 50 - 55 Shore A, tapered to toe 1/4"

431 Blanks:

Top: 3/16" Thermozote
Mid-layer: 1/8" PPT/EVA 30 - 35 Shore A
Base: 1/4" close cell EVA in 50 - 55 Shore A, tapered to toe 1/8"

432 Blanks:

Top: 3/16" Thermozote
Mid-layer: 1/8" PPT/EVA 30 - 35 Shore A
Base: 1/2" close cell EVA in 50 - 55 Shore A, tapered to toe 1/4"

Practitioner Friendly

No Cutting, no gluing, no waiting

Low labor cost

Closed cell material, molds & grinds easily

Ready for vacuum press on mold in 4 - 6 minutes at 350F

Blanks No.	Actual Sizes	Arch Size	Fit Shoe Sizes & Widths
X - Small	L11"x W5"	L4.5"x W2.5"x H3/8"	Sizes 3 - 4.5 / Widths AA - 14E
Small	L12"x W5.5"	L5" x W2.5"x H7/16"	Sizes 5 - 7.5 / Widths AA - 14E
Medium	L13"x W6"	L5.25" x W2.75"x H1/2"	Sizes 8 - 10.5 / Widths AA - 14E
Large	L14"x W6"	L5.7" x W3"x H5/8"	Sizes 11 - 13 / Widths AA - 14E
X - large	L15"x W6.5"	L6" x W3.25"x H11/16"	Sizes 14 - 17 / Widths AA - 14E
XX - Large	L16"x W7"	L6.5" x W3.5"x H3/4"	Sizes 18 - 20 / Widths AA - 14E

Shoe Size/Width Measurement

Because all Apis shoes are made based on American Standard Last measurements, it is easy to figure out the right shoe size, width and depth the patient needs. We have made several ways available for you to follow either method described as below:


1. Electronic measurements: Inputting measurement number at www.apisfootwear.com/findshoessize.asp online tool, and the web site will give you all measurements that is pertinent to the sizes and widths, we highly recommend this tool as it is the most convenient and yields the best fitting results;
2. Using Mt. Emey sizing chart, which you can download from www.apisfootwear.com, detailed instructions are in the chart, or call us to ask free hard copies;
3. Using a standard Ritzstick which Apis supplies at your request;

Using the data base to locate the right shoes. [www.apisfootwear.com / findshoesize.asp](http://www.apisfootwear.com/findshoesize.asp)

Please Log-on Mt. Emey® web-site:

1. Find Shoe Size Web Tool

A. Measure patient's foot data in MM or



B. Input your measurements

Step 2

- Net foot length (heel to longest toe tip)
- Flat ball width (across the ball)
- Ball girth or instep girth
- Click "next"

C. Get the results

Approximately Fitting Shoe Styles, Sizes, and Widths (mm)

Foot Length	Ball Width	Ball Girth	Instep Girth	Heel Girth	Shoe Size	Shoe Width	Available Shoe Styles
260.2	97.5	237.0	259.0	374.0	8.5	2E	9501 9502 9504 9505
260.3	97.5	244.4	258.3	375.0	8.5	2E	502 708 718 888 9912 9901 9921 9951 9925 9923 9928 9931 9938
260.2	98.8	244.8	257.3	363.0	8.5	2E	7021 7022
261.8	98.8	254.8	266.8	385.0	8.5	2E	503 504

* No comparisons have been made to data in black.

2. Using Mt. Emey sizing chart, which you can download from www.apisfootwear.com, detailed instructions are in the chart, or call us to ask free hard copies;

All Mt. Emey & Answer2 shoes/boots are made based on true American Standard Last (FIA Standard). Generic non-branded Ritzstick® or Mt. Emey SizingChart® foot measurement device will be able to tell the size & width. For excessive size & width & depth feet, please download sizing chart from Mt. Emey® technical web-site:

<http://www.bignwideshoes.com/support.html> and fill in the information, & ask us for help. When Printing, use Legal Size Paper, use "do not scale"

3. Using a standard Ritzstick which Apis supplies at your request



Standard Ritzstick measurement

1021

Extremity Equalizer











This device is globally patented for measuring leg-length discrepancy (LLD).
Instructions :

1. barefoot,
2. put books or magazines under the shorter-leg foot
3. place bars on bony landmarks (crests)
4. check horizontal plane for balance
5. check vertical plane for precise balance
6. gradually add lifts to balance leg length
7. measure the thickness of books to decide the lift height
8. write down the height for lift modifications for ordering

Shoe Modifications

Pick any style from the catalog and we will modify the shoes for you ; One low price for each modification.



Apis is the one-stop service station for all of your shoe accommodations and modifications

CUSTOM-MADE ORTHOTICS ORDER FORM

APIS FOOTWEAR 2239 Tyler Ave, S El Monte, CA 91733 Tel: (888)937-2747 Fax: (888)990-2245

P.O. No.: _____

Contact: _____

Bill To: _____ Acct# _____

Address: _____

Ship To: _____

City: _____ State & Zip: _____

Address: _____

City: _____ State & Zip: _____

Ship Via: ☐ UPS 1 day / 2 day / 3 day / Ground ☐ US MAIL ☐ FEDEX 1 day / 2 day / Ground

Patient Name: _____

Shoe Style: _____ Color: _____

Shoe Size (L): _____ (R): _____

Sex: ☐ Male ☐ Female

☐ Check if you want to order shoes

Weight: _____ Height: _____

DIAGNOSIS (diabetic, neuropathy, edema, ankle/

INSERT
SPECIFICATIONS
(quantity, material, etc)



Right Left

Mark prominent areas for off-loading

OTHER SPECIFICATIONS

This form can be downloaded from <http://www.apisfootwear.com/customshoes.htm>

E-Fillable form is available online at <http://apisfootwear.com/default.asp?support.html>

CUSTOM-MADE SHOES ORDER FORM

APIS FOOTWEAR 2239 Tyler Ave, S El Monte, CA 91733 Tel: (888)937-2747 Fax: (888)990-2245

P.O. No.: _____

Contact: _____

Bill To: _____ Acct# _____

Address: _____

Ship To: _____

City: _____ State & Zip: _____

Address: _____

City: _____ State & Zip: _____

Ship Via: ☐ UPS 1 day / 2 day / 3 day / Ground ☐ US MAIL ☐ FEDEX 1 day / 2 day / Ground

Patient Name: _____

Shoe Style: _____ Color: _____

Sex: ☐ Male ☐ Female

Shoe Size (if not custom) (L): _____ (R): _____

Weight: _____ Height: _____


Shoe Closure: ☐ Lace ☐ Velcro

DIAGNOSIS (diabetic, neuropathy, edema, ankle/toe)

LEFT FOOT

Foot length: _____ Ball width: _____

Circumferences: ball _____ instep _____



RIGHT FOOT

Foot length: _____ Ball width: _____

Circumferences: ball _____ instep _____

SHOE MODIFICATIONS

(rocker, lift, velcro/lace, surgical-opening, etc.)

INSERT SPECIFICATIONS

(quantity, material, etc.)



Right Left

Mark areas for off-

loading

OTHER SPECIFICATIONS

Apis Footwear Company

Custom Functional Orthotics Order Form

2239 Tyler Avenue, Unit B, South El Monte, CA 91733 Tel 1-888-937-2747 Fax 1-888-990-2245

P.O No. _____
 Bill To: _____ Account # _____
 Address: _____

Contact: _____
 Ship To: _____
 Address: _____
 City : _____ State _____ Zip _____

City _____ State _____ Zip _____

Phone _____ Fax _____

Ship Via : ☐ UPS ☐ FedEx ☐ ND ☐ 2D ☐ 3D SEL ☐ GRD ☐ US Mail

Patient Name _____ Gender ☐ Male ☐ Female Weight _____ Height _____
 Shoe Style (Refer To Catalog) : _____ Color _____ Shoe Size (Non-custom only) L _____ R _____

CAST MODIFICATION CHART

<p>Modification Methods</p> <p>Modified Root (STD) Inverted _____ Medical Skive _____mm Other _____</p>	<p>Rear Foot Correction</p> <p>Left _____ Degree _____ Right _____ Inverted _____ _____ Everted _____</p>	<p><input type="checkbox"/> Extra Talo-Navic Control <input type="checkbox"/> Extra Mid Arch Control <input type="checkbox"/> Extra Heel Expansion <input type="checkbox"/> Plantar Fascia Groove <input type="checkbox"/> Extra Plaster First Ray <input type="checkbox"/> Cuboid Notch S <input type="checkbox"/> M <input type="checkbox"/> L <input type="checkbox"/></p>
<p>Shell Materials</p> <p><input type="checkbox"/> 4.5 mm Polupro <input type="checkbox"/> 3 mm Polypro <input type="checkbox"/> Carbon Fibre <input type="checkbox"/> EVA - Soft <input type="checkbox"/> EVA - Medium <input type="checkbox"/> EVA - Firm <input type="checkbox"/> Lab Discretion</p>	<p>Shell Selection</p> <p><input type="checkbox"/> Morton's Ext. L: _____mm R: _____mm <input type="checkbox"/> Apertured S M L <input type="checkbox"/> First Ray Cut Out <input type="checkbox"/> Low Bulk Grind Cobra <input type="checkbox"/> Meidal Flare <input type="checkbox"/> Medial Flange Lateral Flange</p>	<p>Posting</p> <p><input type="checkbox"/> Rear Foot EVA posts <input type="checkbox"/> Heel Raise L: _____mm R: _____mm <input type="checkbox"/> No Heel Plate <input type="checkbox"/> EVA Arch Fill- Density _____ <input type="checkbox"/> Forefoot EVA Posts -Left _____Degree Inv/Eversion -Right _____Degree Inv/Eversion</p>
<p>Cover Length</p> <p><input type="checkbox"/> Shell <input type="checkbox"/> Web <input type="checkbox"/> Full</p> <p>Met Domes</p> <p><input type="checkbox"/> Small <input type="checkbox"/> Medium <input type="checkbox"/> Large <input type="checkbox"/> Extra Large <input type="checkbox"/> Others _____</p>	<p>Padding & Cover</p> <p>Cover Materials</p> <p><input type="checkbox"/> Vinyl (STD.) <input type="checkbox"/> Leather <input type="checkbox"/> Neoprene <input type="checkbox"/> Other _____</p>	<p>Padding</p> <p>Shell Forefoot</p> <p><input type="checkbox"/> None <input type="checkbox"/> <input type="checkbox"/> 1.5 mm PPT <input type="checkbox"/> <input type="checkbox"/> 3.0mm PPT <input type="checkbox"/></p>

Notes :

Mark Offload Area



Time _____

Apis Footwear Order Form

Toll Free Tel: 1-888-937-2747 (YES-APIS); Toll Free Fax: 1-888-990-2245 E-mail: orders@apisfootwear.com

Order Date: _____ : Order Taken by: _____ Input-in. by: _____ **Order NO:** _____

Sold to: **Account No:** _____

Authorized Ordering Person : _____ (First name) _____ (Last Name)

E-mail Address to E-Notify the E-Billing status: _____

Telephone No: (_____) _____ Fax No.: (_____) _____

Bill to **Account Name:** _____

Add: _____, City: _____ State _____ Zip _____

Ship To: (Same as Bill to) (Different from Bill To)

Ship to Name: _____

Add: _____, City: _____

Telephone No: (_____) _____ State _____ Zip _____

Customer's Reference PO No. : _____

Ship By: UPS _____ FEDEX ___ U.S MAIL ___ DELIVER ___ PICK UP ___ OTHER _____

Item #	Color	Width	Size	QTY	Left Size	Right Size	Patient's Name

Please Process As: _____ (Regular Daily Ship) ; _____ (Priority Ship Today); _____ (Extremely Urgent, Must Ship Today)



Apis Footwear Return Policy

90-Day Satisfaction Guarantee

Your satisfaction is our goal! We strive to provide our customers with the highest level of service possible. **Before returning the product for refund, please make sure:**

- (1) All products must be in their original, unworn condition. Please try on the product on carpeting to prevent scuffing.
- (2) All products must be returned in their original package and returns must be placed inside a packing box to avoid damage to the product and packaging.
- (3) All accessories (if any) came with the product should be returned in their original, unworn condition.
- (4) All returns must be made within 90 days of the issuance of invoice.
- (5) Please complete the return items section below and include it with your return order.
- (6) Shipping costs for return of damaged goods and subsequent replacement goods will be paid Commerce Atlantic Corporation dba Apis Footwear.
- (7) If a shipping error occurs due to Commerce Atlantic Corporation dba Apis Footwear, the company will pay freight both ways for return/exchange of goods. If an ordering facility error occurs, the ordering facility will pay freight both ways without a restocking charge.

There is no need to contact our customer service for a return unless you have a question or need assistance. Please save your tracking number.

Return For Refund

Refunds for returned product will be issued within one week of receipt of the item/items. After your credit is processed, it may take an additional week for the funds to show in your account. Most banks take approximately a week to transfer funds on a credit card.

Returned Items

Return Reason	Item Number & Description	Return Quantity

Question? Call Toll-free **1-888-937-2747** (8am to 5pm PST Mon.-Fri.)

Apis Footwear Company, 2239 Tyler Avenue, Unit B, South El Monte, CA 91733
www.apisfootwear.com email: sales@apisfootwear.com

*9708 Walking Shoes
D, 4E, 6E & 9E,
7, 7.5, 8, 8.5, 9, 9.5, 10, 10.5, 11,
11.5, 12, 13, 14, 15, 16 & 17*



Apis Footwear Company

2239 Tyler Avenue, South El Monte, CA 91733 (West Coast Distribution Center)

13026 Forest Center Court, Louisville, KY 40223 Tel: 502-245-2747 (East Coast Distribution Center)

Tel: 888-937-2747 (CA) Tel: 1-888-777-0448 (KY) Fax: 888-990-2245

www.apisfootwear.com

Apis
Footwear Co.