



Orlando® KAFO Pre-assembled Femoral Fx Orthosis Pre-fabricated and pre-assembled for many full leg orthopaedic indications. Especially useful for distal femoral or proximal tibial fractures or hemiplegic patients following CVA. Narrow M-L, longer A-P thigh design for more anatomical fit. Free motion, adjustable, or drop-lock knee joints allow for specific range-of-motion and functional knee control. Single axis or solid ankle shoe insert with extended footplate will prevent flaccid foot drop for post-CVA patients. Plastic joint covers for cosmesis and support. Components also available separately. Available with pelvic component, please see HKAFO on page 44.

- Trimmable polyethylene components offer more comfortable total contact
- Pre-assembled—not a kit
- Lightweight, durable polyethylene femoral and tibial components with a low density tongue
- Excellent for office or clinical application
- Available with adjustable joints

	Proximal Circ.	Mid Thigh Circ.	Distal Circ.	Largest Calf Circ.	Length (Perineum to floor)	U.S. Shoe Size
Small	17"-21"	14"-19"	13"-16"	11"-12"	28"-31"	Men 5-7 Women 6 or less
Medium	20"-24"	19"-24"	15"-18"	13"-15"	31"-34"	Men 7½-10 Women 6½ & up
Large	23"-27"	23"-27"	17"-20"	15"-17"	34"-37"	Men 10½-13

All components interchangeable as well as available separately. Refer to page 49 for component sizing and details.

Free Motion KAFO With 1/8" free motion knee joints, free ankle motion shoe insert

	Left Model #	Right Model #
Small	1690	1695
Medium	1691	1696
Large	1692	1697

Adjustable Range-of-Motion KAFO

With 1/8" adjustable range-of-motion knee joints, free ankle motion shoe insert

	Left Model #	Right Model #
Small	1690.01	1695.01
Medium	1691.01	1696.01
Large	1692.01	1697.01

Drop Lock KAFO For hemiplegic indications with 1/8" drop lock knee joints, solid ankle shoe insert

	Left Model #	Right Model #
Small	1700	1705
Medium	1701	1706
Large	1702	1707