

CLOSED-END PVA SIZING CHART

Closed-End PVA Size	Proximal Measurement	Distal Measurement
4 - 0	6" - 9"	2"
6 - 1	9" - 12"	3"
6 - 3		4" - 6"
6 - 5		7" - 10"
8 - 3	12" - 16"	4" - 6"
8 - 5		7" - 10"
10 - 3	16" - 20"	4" - 6"
10 - 5		7" - 10"
10 - 7		11" - 14"
12 - 3	20" - 24"	4" - 6"
12 - 5		7" - 10"
12 - 7		11" - 14"
14 - 3	24" - 28"	4" - 6"
14 - 5		7" - 10"
14 - 7		11" - 14"
16 - 5	28" - 32"	7" - 10"
16 - 7		11" - 14"

EXAMPLE: If your proximal size is 15", and the distal size is 9", then the correct Closed-End PVA size would be an 8-5.