



The Miami JTO® Thoracic Extension is designed to be used exclusively with any Miami J® collar to create a CTO for treating low cervical-high thoracic injuries or for more control than a collar alone. The floating back helps maintain desired alignment when the patient is supine. The modular chin strut enhances flexion control by directly bridging the chin and sternum. Free of metal for imaging compatibility. Breathable material wicks moisture away. Easy Velcro® attachment.

Miami JTO® Thoracic Extension

Order No.	Description
MT-900	Miami JTO® with Pads
MTR-900	Miami JTO® with Extra Pads
MTP-950	Replacement Pads
MT-900 FRT	Front Only
MT-900 BK	Back Only