

UltraSeal Liner (Low Ring) Sizing Chart

Ring: The low ring begins 5 inches from the center of the distal end.

How to Size: Measure the circumference of the residual limb at 6 cm above the distal end.

MEASURED SIZE (CM)	LOW RING LOCKING 3.0MM UNIFORM	LOW RING CUSHION 3.0MM UNIFORM
20-21 cm	USLSL20	USLSC20
22-23 cm	USLSL22	USLSC20
24-25 cm	USLSL24	USLSC24
26-27 cm	USLSL26	USLSC26
28-29 cm	USLSL28	USLSC28
30-31 cm	USLSL30	USLSC30
32-34 cm	USLSL32	USLSC32
35-37 cm	USLSL35	USLSC35
38-40 cm	USLSL38	USLSC38
41-43 cm	USLSL41	USLSC41
44-47 cm	USLSL44	USLSC44