

3 Materials

3.1 Plaster Fabrication

AK Casting Brims

Version “CASCD”

With fixation holes for AK Brim Stand

- More rigid than Version “Plaster Model”
- Sold as single unit or as complete set
- 4 different socket model available

Size Recommendation

Measure the circumference under the Tuber os ischii.
Choose the brim in accordance to the stump condition
(i.e. muscles) one to two sizes smaller.



MM220 Soft (quadrilateral) Per single unit

Order No.	Article No.
MM220-(Size)(Side)	4 304 070 0X 00 0ZZ
Available sizes in left and right (in cm)	42, 44, 46, 48, 50, 52, 54, 56

X = Side: 1 = Left, 2 = Right _ ZZ = Size (check available sizes)



Complete set (all sizes, left and right)

Order No.	Content	Article No.
MM220	16 Units	4 304 070 03 00 000

MM221 Hard (quadrilateral) Per single unit

Order No.	Article No.
MM221-(Size)(Side)	4 304 080 0X 00 0ZZ
Available sizes in left and right (in cm)	28, 30, 32, 34, 36, 38, 40, 42, 44, 46, 48, 50, 52, 54, 56

X = Side: 1 = Left, 2 = Right _ ZZ = Size (check available sizes)



Complete set (all sizes, left and right)

Order No.	Content	Article No.
MM221	30 Units	4 304 080 03 00 000

3.1 Plaster Fabrication

MM222 CNC (quadrilateral) Per single unit

Order No.	Article No.
MM222-(Size)(Side)	4 304 090 0X 00 0ZZ
Available sizes in left and right (in cm)	36, 38, 40, 42, 44, 46, 48, 50, 52, 54, 56, 58, 60, 62, 64

X=Side: 1=Left, 2=Right _ ZZ=Size (check available sizes)

Complete set (all sizes, left and right)

Order No.	Content	Article No.
MM222	30 Units	4 304 090 03 00 000

**MM223 Ischial Containment (modified quad/IC) Per single unit**

Order No.	Article No.
MM223-(Size)(Side)	4 304 100 0X 00 0ZZ
Available sizes in left and right (in cm)	36, 37, 38, 39, 40, 41, 42, 43, 44, 45, 46, 47, 48, 49, 50, 52, 54, 56, 58, 60, 62, 64

X=Side: 1=Left, 2=Right _ ZZ=Size (check available sizes)

Complete set (all sizes, left and right)

Order No.	Content	Article No.
MM223	44 Units	4 304 100 03 00 000

**AK Portable Brim Stand**

- As a fitting stand for transfemoral trial or final sockets
- For immediate postoperative use
- Simple and individual length adjustable
- Easy attachment of the AK Casting Brims

Order No.	Article No.
MM200-1	4 302 010 00 00 000

