

Comfor™ Clavicle Brace - Hook and Loop Closure

Application Instructions



Comfor™ Clavicle Brace - Hook and Loop Closure

- Our most popular hook and loop closure clavicle brace
- Stockinette covered felt
- Comfortable wide padding for proper immobilization
- Felt back pad helps prevent irritation
- Hook and loop closure for easy application

1. Remove brace from package and loosen hook and loop part of strap.
2. Feed strap through plastic buckle opening, adjust according and tighten by attaching hook and loop strap.

Reorder No.	Size	Description
0814 1230	XX-Small	Infant
0814 1231	X-Small	Child
0814 1232	Small	Child/Youth
0814 1233	Medium	Youth/Small Adult
0814 1234	Large	Adult
0814 1235	X-Large	Large Adult

BIRD & CRONIN INC.

Designers and Manufacturers of Quality Medical Products

1200 Trapp Road • Eagan, Minnesota USA 55121

651/683-1111 • Fax 651/688-9855 • 800/328-1095 • Fax 800/279-7934

e-mail sales@birdcronin.com • www.birdcronin.com

©2005 Bird & Cronin, Inc. All Rights Reserved - 10001230