

### WHAT'S IN THE BOX

Alpha Hybrid AK Sleeve  
Prosthetist Instructions  
Patient Instructions

Please provide the amputee with the **Patient Instructions** included in the box.

### CUTTING

It is not likely that an Alpha Hybrid AK Sleeve will have to be cut to a shorter length. If trimming is required, however, use large shears, a paper cutter, or a Hand-Held Cutting Wheel (Part No. 700-TL002). If using the Hand-Held Cutting Wheel, place the sleeve on a cutting board or a thick piece of plastic. Using a straight edge as a guide, roll the cutting wheel across the sleeve. Cutting the sleeve too short is not covered under warranty.

### WARRANTY

WillowWood warrants each Alpha Hybrid AK Sleeve from the date of invoice for a period of ninety (90) days against defects in material and workmanship.

A product returned during the warranty period will be evaluated and repaired or replaced at WillowWood's discretion. To avoid restocking fees, you must obtain a Return Authorization (RA) Number prior to the return of any WillowWood product. Call 740-869-3377 or 800-848-4930 to obtain an RA.

There is no warranty, expressed or implied, for damage caused by abuse, mishandling or accident. The warranty is void if the fitting and routine care instructions are not followed. For your protection and ours, please thoroughly clean all used products prior to return.

### Warranty Disclaimer

WillowWood warrants that each product manufactured will, at the time of delivery, be of workmanlike quality and substantially free of defects. WILLOWWOOD MAKES NO OTHER WARRANTY, IMPLIED, OR EXPRESSED, AND MAKES NO WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. This warranty shall terminate immediately upon an action to combine our products with other materials or in any manner to change the nature of our products. The sole remedy is replacement of the products or credit for the products. Willow Wood's liability shall not exceed the purchase price of the product. WillowWood shall not be liable for any indirect, incidental, or consequential damage.

### WillowWood Retention of Rights

WillowWood retains all intellectual property rights reflected or incorporated in its physical products, regardless of the transfer of the physical products to another party or parties.



**Alpha Products can melt or burn if exposed to high temperatures or flame. Do not expose your Alpha product to these conditions.**



**This Alpha product is intended for use on a single patient. Use of the product with multiple patients could lead to cross contamination between patients.**



**WARNING: To avoid danger of suffocation, keep this product away from babies or children.**

# WillowWood<sup>®</sup>

The Ohio Willow Wood Company  
15441 Scioto Darby Road  
Mt. Sterling, OH 43143  
phone 740.869.3377 / 800.848.4930  
fax 740.869.4374 [www.willowwoodco.com](http://www.willowwoodco.com)



Ohio Willow Wood Company B.V.  
Keizersgracht 62/64  
1015 CS Amsterdam  
The Netherlands

