



ottobock.

Skeo 3D

For freedom of movement and protection

Quality for life



Discover new wearer comfort!


Haven't you always wanted a liner that is based on the anatomy of your lower leg and thus makes wearing it more comfortable?


This is what a liner with the Anatomy-Specific Geometry (ASG for short) offers you. The material is thicker in sensitive areas to provide extra protection, and thinner in other areas for greater flexibility and support.


6Y77 Skeo 3D liner –


Your benefits:


-  **Thin material**
-  **Thick material**


 **Patella/Kneecap**
Less pressure on the patella

 **Surrounds the patella**
Stabilizing support

 **Fibular head**
Protects the sensitive region

 **Shinbone/Tibia**
Protection from impacts

 **Distal cushion**
Protects the sensitive residual limb end



Anatomy-Specific Geometry (ASG)

The liner features different wall thicknesses. The material is thicker in sensitive areas to provide protection and thinner in other areas for greater flexibility.

Special inner coating

A silky inner coating ensures greater comfort without sticking.

Increased knee flexion (20°)

The liner is pre-flexed so the knee bends easily. As there is no excess material to create wrinkles in the hollow of the knee, uncomfortable pressure points are avoided.


Donning Aid

The logo is useful as a donning aid for a better orientation for the patient.

Distal connector

Soft distal cap with connector



 Easier to put on and take off



 MagnoFlex Lock with flexible pin

Our fitting recommendation for you:

Due to their good suspension properties, silicone liners such as the 6Y77 Skeo 3D are recommended for use in combination with a closure system such as the 6A40 MagnoFlex Lock. The MagnoFlex Lock has a special flexible pin and the innovative shuttle lock with integrated magnet making to guide the pin easier than ever before to join the liner and prosthesis. The ideal lower limb fitting for residual limbs with good soft tissue coverage.