

## REBOUND® POST-OP BRACING EASY TO FIT. EASY TO USE.

The Rebound Post-Op Knee & Rebound Post-Op Elbow are designed with two priorities in mind: ease of fitting for practitioners and ease of use for patients.

## Simple and intuitive.

Tool-free design makes fitting fast and easy, while the numbered strapping order and color-coded patient touch points simplify patient education and product use. Just tell your patients, "Blue is You™."





## Total control.

Use the OneTouch™ hinge to quickly and easily set patient range-of-motion. Then, secure your range-of-motion settings in place with lockout clips to further reduce the risk of patient non-compliance.

### Comfortable and secure.

Lighter than the competition, Rebound Post-op Knee and Rebound Post-Op Elbow have ergonomic paddles and contoured hinge struts that improve comfort and anatomical fit, as well as Anti-Migration Features to help ensure secure suspension.



## Something for everyone.

The Rebound Post-Op bracing series ensures you have everything you need to fit any patient you encounter. Model options include:

- Universal versions of both the knee and elbow are suitable for most users
- For the knee; Pediatric, XL or left/right Contoured model accommodates when a more anatomical fit is required



Rebound® Post-Op Knee Contoured



Knee Universal & XL



Knee Pediatric



Rebound®Post-Op Rebound®Post-Op Rebound®Post-Op Elbow

# REBOUND® POST-OP ELBOW



### Features and Benefits

- OneTouch™ Drop-lock hinge
- 2 Four telescoping paddles for easy adaptability to individual post-op protocols
- 3 Strap guides keep straps in place
- 4 Moldable stability cuffs promote user comfort
- 5 Blue is You<sup>™</sup>- color-coded buckles and drop lock for easy patient education and product use

## REBOUND® POST-OP KNEE

### Features and Benefits

- Large buckles feature hook closures and grip-friendly ridges for simple brace application and removal
- 2 Numbered frame for patient ease-of-use
- 3 Blue is You<sup>™</sup> Color-coded drop lock and buckles for easy patient education and product use
- 4 Lock-out ROM pull tabs offer 10° to 120° of ROM control, with OneTouch™ drop-lock from 10° to 30°
- 5 Anti-Migration System wrap for secure suspension
- 6 Quick-fit buckle with a lever arm helps ensure optimal strap tension and brace suspension
- 7 Telescoping frame features easy-to-engage OneTouch™ button with secure locking feature
- 8 Alligator strap tabs enable easy strap trimming
- Moldable, anatomically contoured struts help ensure a precise fit



#### **CONTINUUM OF CARE**

Össur offers a breadth of Injury Solutions that follow a recommended continuum of care model. After knee injury or surgery, protect and treat by restricting range-of-motion with the Rebound® Post-Op Knee, while managing pain with Össur cold therapy solutions. In some cases, rehab with a Functional Healing solution may be necessary, such as with the Rebound® ACL or Rebound® Cartilage. Finally, return to pre-injury activity with CTi® Mission or Rebound® DUAL. Choose Össur products as rehabilitation progresses to follow the Össur Continuum of Care.

### POST-OP



Rebound® Post-Op Knee Controlled range-of-motion in a light-weight, quick-fitting design.



Cold Rush® Cold Therapy Convenient, whisper quiet and easy-to-use, Cold Rush delivers hours of cold therapy when and where you need it.

## **REHAB**



Rebound® PCL & Rebound® ACL Designed to optimize rehab by reducing load on the PCL/ACL respectively, by applying a dynamic, anatomicallycorrect force during nonsurgical or postsurgical treatment.

#### Rebound® Cartilage A protective, functional solution designed to support regeneration of the knee cartilage after injury or surgery.

## 3 RECOVERY



CTi® Mission Comfortable, lightweight and durable ligament knee brace, providing stabilization of the knee joint and versatile for both everyday activities and sports.



Rebound® DUAL Functional support for ligament instabilities, including ligament instability with osteoarthritis, in a single, versatile aluminum brace. Also available in step-through design.

## ORDERING INFORMATION

#### **REBOUND POST-OP KNEE**

Size	Thigh Circumference*	Telescoping	Cool	Full Foam
Pediatric	12"-25" (30.5cm - 63.5cm)	15.5"-21.5" (39.5 - 54.5cm)	B-8381470	N/A
Universal	16"-32" (41cm - 81cm)	18.5"-27.5" (47 - 69.5cm)	B-8381480	B-8382480
Left contoured	16"-32" (41cm - 81cm)	18.5"-27.5" (47 - 69.5cm)	B-8381680	B-8382680
Right contoured	16"-32" (41cm - 81cm)	18.5"-27.5" (47 - 69.5cm)	B-8381580	B-8382580
Universal XL	16"-38" (41cm - 96.5cm)	18.5"-27.5" (47 - 69.5cm)	B-8381490	B-8382490

<sup>\*</sup>Measurement taken at 6" (15cm) above mid-patella

#### REPLACEMENT PADS

Size	Cool	Full Foam	AMS Accessory Kit		
Pediatric	B-8381072	N/A	B-8381093		
Universal†	B-8381082	B-8382082	B-8381094		
Universal XL	B-8381092	B-8382092	B-8381095		

<sup>†</sup>For universal, left contoured, right contoured

#### KNEE ACCESSORY KITS

Accessories	Part Number
Knee ROM Stop locks Kit	B-8380005
Incision Pad Kit (Pack of 10)	B-703180005+
Standard pad and buckle kit	B-8381096

Visit ossur.com/rebound for more information.

#### **REBOUND POST-OP ELBOW**

Side <sup>‡</sup>	Part Number
Left	B-8391680
Right	B-8391580

±Universal sizing

#### FLBOW ACCESSORY KITS

Description	Part Number	
Elbow Liner Kit	B-8390002	
Elbow Wrist Control Kit	B-8390003	
Elbow Sling Kit	B-8390004	
Elbow ROM Stop Locks Kit	B-8390005	
Strap Tab Kit	B-8390006	
Adhesive Hook Sticker Kit	B-8390007	
Anti-Migration Pad Kit	B-8390008	
Elbow Strap Kit	B-8390009	

FOLLOW ÖSSUR ON











USA (800) 233-6263 WWW.OSSUR.COM CANADA (800) 663-5982 WWW.OSSUR.CA



