

Chip Pad Dorsum

Ideal for use on dorsum of hand or foot. Helps with softening scars from carpal tunnel surgery, or underneath compression bandages over pockets of lymphedema. Fits into Dorsal Pocket Glove.

Size	Single
Small	2231-DO (8cm x 9cm)
Large	2233-DO (9cm x 14cm)



Chip Pad Axilla

Use under arm for treatment of lymphedema, especially over areas of post-operative scarring.

Size	Single
Small	2231-AX (10cm x 19cm)
Large	2233-AX (18cm x 27cm)



Chip Pad Radial Half

Radial pathways create low pressure channels to promote fluid drainage. Use over effected areas or pockets of edema.

Size	Single
One size fits all	2230-RH (15cm x 25cm)



Chip Pad Radial Full

Versatile, effective design is ideal for post-surgical scarring and/or pockets of lymphedema.

Size	Single
Small	2231-RF (18cm diameter)
Large	2233-RF (23cm diameter)



Chip Pad Bra

Ideal for post mastectomy/lumpectomy edema. Pair with Comprevest.

Size	Single
Black	
Half	2230-BRH
Full	2230-BRF
White	
Half	2280-BRH
Full	2280-BRF



Chip Pad Bra Half




Chip Pad Bra Full

Chip Pads Genifit™

The Genifit series are gender specific Chip Pads, designed for the treatment of genital edema. Each Chip Pad includes open-cell foam chips sewn in channel formations to assist in massaging tissue and facilitate drainage. The Chip Pads are made with moisture-wicking fabric and may be machine washed and dried for repeated use.

Male



Tall

The rise of the tall chip pad ranges in length from 32–44cm. See the measurement guide for reference.




Short

The rise of the short chip pad ranges in length from 14–22cm. See the measurement guide for reference.

Male Genifit					
	Small	Medium	Large	X Large	XX Large
Width	28cm	32cm	36cm	40cm	44cm
Scrotal Circ.	26cm	32cm	38cm	44cm	51cm
Rise – Short	14cm	16cm	18cm	20cm	22cm
Rise – Tall	32cm	35cm	38cm	41cm	44cm
Short	2231-GMS	2232-GMS	2233-GMS	2234-GMS	2235-GMS
Tall	2231-GMT	2232-GMT	2233-GMT	2234-GMT	2235-GMT

Female



Tall

The rise of the tall Genifit female chip pad ranges in length from 32–41cm. See the measurement guide for reference.



Short

The rise of the short female chip pad ranges in length from 14–20cm. See the measurement guide for reference.

Female Genifit				
	Small	Medium	Large	X Large
Width	28cm	32cm	36cm	40cm
Rise – Short	14cm	16cm	18cm	20cm
Rise – Tall	32cm	35cm	38cm	41cm
Short	2231-GFS	2232-GFS	2233-GFS	2234-GFS
Tall	2231-GFT	2232-GFT	2233-GFT	2234-GFT