



ContexGel.soft

ContexGel-Liner, anterior wall thickness 6 mm, posterior wall thickness 3 mm, made of soft ContexGel, with outer matrix

Designated use

- activity level (walky): 1 to 3
- for below-knee amputees
- suitable to substitute soft tissue at the tibia edge and stump end
- for bony and difficult leg stumps

Characteristics

- thermoplastic mouldable at 60 °C (140 °F)
- optimal wearing comfort and very good shock absorption in the anterior stump area
- especially soft functional textile cover for easily rolling the liner up or down
- outer matrix limits longitudinal stroke in the distal area
- thinner wall thickness in the posterior dorsal area improves knee flexion at prosthetic socket
- functional textile cover minimises longitudinal expansion
- excellent circular expansion adjusts to volume changes at the residual limb
- distal rotation control supports guidance of prosthesis
- available with distal connections for all Streifeneder pin types and locking systems
- available without distal connection for vacuum system technology (VST)

Material


- anterior wall thickness 6 mm, posterior wall thickness 3 mm, conical shape
- soft ContexGel with medical white oil, antibacterially fitted
- approx. 25 – 30 Shore 00
- abrasion-proof textile cover made of polyamide/elastane
- item 26L100 with PUR-rotation control
- colour: anthracite

Warranty

- 6 months



Dimension sketch:



* 4,5 mm without distal connection
14 mm with distal connection

total length: approx. 420 mm

Size	Distal	Proximal	Distal connection	Item-No.
S	16 - 23 cm	25 - 36 cm	with	26L100/S
M	19 - 30 cm	30 - 43 cm	with	26L100/M
L	24 - 34 cm	38 - 58 cm	with	26L100/L
XL	28 - 39 cm	41 - 63 cm	with	26L100/XL
XXL	32 - 44 cm	43 - 69 cm	with	26L100/XXL
S	16 - 23 cm	25 - 36 cm	without	26L200/S
M	19 - 30 cm	30 - 43 cm	without	26L200/M
L	24 - 34 cm	38 - 58 cm	without	26L200/L
XL	28 - 39 cm	41 - 63 cm	without	26L200/XL
XXL	32 - 44 cm	43 - 69 cm	without	26L200/XXL

Scope of supply does include stand.

