



Extender Plus and Extender Post-Op Knee Braces



Telescopes from 22" – 30" (56 – 76 cm) in length and accommodates circumference up to 34" (86 cm). Available options include cuffs or a short version which adjusts from 18" – 26" (46 – 66 cm) in length.

Common Examples of Use

- Post-op rehabilitation
- Soft tissue injuries of the knee
- ACL, PCL, MCL, LCL repairs / injuries
- Tibial plateau fractures
- Osteochondral repairs
- Meniscal repairs
- Patella tendon repairs
- Condylar fractures
- Sprains / strains of the knee
- HTOs (High Tibial Osteotomy)
- Stable or internally fixed fractures of the knee
- Shoe insert is indicated when support or immobilization of the foot/ankle is needed

Product Features

- Versatile and quick application with push button length adjustment
- Ankle plate flared out for comfort
- Quick release buckles
- Optional shoe insert

Category: Post Op.

Instructions for Use

Instructions for Use

Part Numbers

Part Numbers

Part #	Description
EK087000	Extender Plus
EK088000	Extender Plus, w/ Cuffs
EK009000	Extender Plus, Short
EK0071XX	Extender Plus Shoe Insert, Left
EK0072XX	Extender Plus Shoe Insert, Right
SK587000	Extender w/ Quicklock
SK587010	Extender w/ Shear Force Straps and Quicklock
SM0071XX	Extender Shoe Insert, Left
SM0072XX	Extender Shoe Insert, Right

Product Sizing Chart

Shoe Insert Sizing Chart				
Size #	Size	Men's	Women's	Euro
XX = 03	S	Up to 4	Up to 5	Up to 34.5
XX = 05	M	5 – 9	6 – 10	35 – 41
XX = 09	L	10 – 14	11 – 15	41.5 – 47.5