# FLEX-RUN™





The Flex-Run with Nike Sole offers unmatched performance and efficiency for amputees and prosthetists. Össur and Nike have partnered to design a running foot from the ground up, leveraging Össur's expertise in running prosthetics with Nike's world renowned traction and sole technologies.

The improved Flex-Run features a longer toe lever and efficient energy return. The Nike Sole offers great traction to help improve performance for the runner.

Flex-Run with Nike Sole is appropriate for individuals OR users with transfemoral and transtibial amputation participating in high impact activities such as recreational jogging, trail running, distance running and triathlons. It is available in both Standard and Junior (LP) models.



User Information					
Amputation Level:	Transfemoral and Transtibial				
Impact Level:	High to Extreme				
Maximum Patient Weight:	130kg (287lbs)				



$  \triangle_{\wedge}  $			
0			
WATERPROC	)F		

Foot information					
Categories:	1-9				
Sizes:	One Size Per Category				
Weight of Foot (Cat. 5):	400g (14.1oz) 694g (24.4oz) with Sole and Male Pyramid				
Build Height:	261mm (10 1/4") 277mm (10 7/8") with Sole				
Adapter Options:	Male or Female Bolt on Pyramid				

# ORDERING INFORMATION

#### FOR COMPLETE FLEX-RUN KIT WITH MALE PYRAMID



#### FOR COMPLETE FLEX-RUN KIT WITH FEMALE PYRAMID






# FLEX-RUN™

## CATEGORY SELECTION CHART

Weight kg	37-44	45-52	53-59	60-68	69-77	78-88	89-100	101-116	117-130
Weight lbs	81-96	97-115	116-130	131-150	151-170	171-194	195-220	221-256	257-287
High to Extreme Impact Level	1	2	3	4	5	6	7	8	9


## **CLEARANCE**



Note: Sole adds 10 mm (3/8") to build height.

### KITS


#### FLEX-RUN KIT



Note: Flex-Run is excluded from trial period. For LP Flex-Run see Flex-Run Junior. For information on special or custom orders, please contact your local Össur representative.

#### Includes:

- Foot Module
- Attachment Plate (can be supplied bonded on the foot upon request)
- Nike Sole
- Bag for Nike Sole



# FLEX-RUN™

#### KITS

#### NIKE SOLE KIT

Part#	Size	Cat.
CM150163	One Size	1-9

#### Includes:

- Nike SoleMedallion
- LoctiteBag for Nike Sole



#### FOOT ADAPTER SETS

Picture	Part#	Description	Weight	Bolt Size	Cat.
1.	FWX00030	Male Pyramid Foot Adapter Set	75g	M10 x 40mm	1-6
2.	CM150161	Female Pyramid Foot Adapter Set	100g	M10 x 35mm	1-6
2.	CM150162	Female Pyramid Foot Adapter Set	186g	M12 x 50mm	7-9





#### REPLACEMENT BOLT

Part#	Size	Cat.
CM130221	M10 x 40	1-6

# **COMPONENTS**

Part#	Description	Cat.
FWX40001	Attachment Plate	1-6
FWX40002	Attachment Plate	7-9