


How to measure your knee for a **GenuTrain® S Pro** knee support.



Measure the circumference of your leg at two points:

4 $\frac{3}{4}$ " below the knee and 5 $\frac{1}{2}$ " above the middle of the knee.

Size	A	B
	Circumference (inch)	Circumference (inch)
1	11-12 $\frac{1}{4}$	15-16 $\frac{1}{4}$
2	12 $\frac{1}{4}$ -13 $\frac{1}{2}$	16 $\frac{1}{4}$ -17 $\frac{1}{4}$
3	13 $\frac{1}{2}$ -14 $\frac{1}{2}$	17 $\frac{1}{4}$ -18 $\frac{1}{2}$
4	14 $\frac{1}{2}$ -15 $\frac{3}{4}$	18 $\frac{1}{2}$ -19 $\frac{3}{4}$
5	15 $\frac{3}{4}$ -17	19 $\frac{3}{4}$ -20 $\frac{3}{4}$
6	17-18	20 $\frac{3}{4}$ -22
7	18-19 $\frac{1}{4}$	22-23 $\frac{1}{4}$