

## STRUCTURE



1. Three-dimensional Train active knit - for the ideal fit and optimal neuromuscular joint stabilization (proprioception)
2. Adjustable straps - allows for an individualized level of pressure for better fit and therapy
3. Joint side bars - Straps connect directly to side bars for superior stability of the knee
4. Viscoelastic pad surrounds the patella and facilitates correct positioning and pain reduction
5. Soft knit fabric - Breathable knit is comfortable to wear and keeps moisture away from the body
6. Reduced pressure at the edges - prevents constriction