



# SOCKET ATTACHMENT DEVICE



1K01



1K02

## • Wooden Blocks

These socket attachment wooden blocks are bonded at the distal extremity of wood sockets. They can be assembled to any screwing system with European 4-hole pattern.

REF.	Description	Max. Patient's Weight	Material	Weight	TH	Circumference
1K01	Medium Model	100 kg	Wood	240 g	53 mm	128 x 128 mm
1K02	Large Model for Trans-Femoral Socket with European 4-Hole Pattern (Off Center)	100 kg	Wood	410 g	53 mm	180 x 180 mm



## • Polyurethane Linking Block for Socket with European 4-Hole Pattern.

1K25 linking block is glued at the distal end of endoskeletal laminated sockets.

REF.	Max. Patient's Weight	Material	Weight	TH	BH
1K25	100 kg	Polyurethane	280 g	36 mm	36 mm



## • Junction Cup for Lamination or Thermoforming

For lamination or thermoforming. This aluminium cup allows connection of any component with european 4-hole 36 mm pattern.

REF.	Max. Patient's Weight	Material	Weight	TH	BH
1K90	125 kg	Aluminium	50 g	22 mm	6 mm



## • Child Junction Cup for Lamination or Thermoforming

For lamination or thermoforming. This aluminium cup allows to connect components with 4-hole 30 mm pattern to the distal extremity of the child socket.

REF.	Max. Patient's Weight	Material	Weight	TH	BH
1K95	45 kg	Aluminium	45 g	22 mm	5 mm

