

SOCKET ATTACHMENT DEVICE



• Wood Attachment Block

This socket attachment wooden block is bonded at the distal end of wooden sockets.

It can be assembled to any screwing system with European 4-hole pattern.

Item N°.	Description	Max. Patient's Weight	Material	Weight	TH	Circumference
1K01	Medium Model	100 kg	Wood	240 g	53 mm	128 x 128 mm



• Polyurethane Linking Block for Socket with European 4-Hole Pattern

1K25 linking block is glued at the distal end of endoskeletal laminated sockets

Item N°.	Max. Patient's Weight	Material	Weight	TH	BH
1K25	100 kg	Polyurethane	280 g	36 mm	36 mm



• Junction Cup with European 4-Hole Pattern for Lamination or Thermoforming

This aluminium cup allows to connect components to the distal end of the socket.

Item N°.	Max. Patient's Weight	Material	Weight	TH	BH
1K90	125 kg	Aluminium	50 g	22 mm	6 mm
1S108 Adapter for 1K90	-	Plastics	15 g	17,5 mm	17,5 mm



• Child Junction Cup for Lamination or Thermoforming

For lamination of thermoforming. This aluminium cup allows to connect components with 4-hole 30 mm pattern to the distal end of the child socket

Item N°.	Max. Patient's Weight	Material	Weight	TH	BH
1K95	45 kg	Aluminium	45 g	22 mm	5 mm