


circaid

reduction kit sizing information for leg



circaid reduction kit - upper leg			
Circumference in cm	regular	wide	
cG	< 80	< 125	
cE2	< 65	< 110	
Length in cm	short	standard	long
∠ D1 - G	< 37	37 -41	> 41

Each upper leg kit includes:

- 1 reduction upper leg component
- 1 reduction shelf strap
- 6 hook and loop stays
- 1 Built-In-Tension system guide card
- 1 paper measuring tape
- 1 direction for use

circaid reduction kit - lower leg			
Circumference in cm	regular	wide	
cC	< 65	< 80	
cB	< 55	< 65	
Length in cm	short	standard	long
∠ mall - D1	< 32	32-36	> 36

Each lower leg kit includes:

- 1 reduction lower leg component
- 1 pair reduction kit undersleeve leg
- 1 reduction shelf strap
- 6 hook and loop stays
- 1 Built-In-Tension system guide card
- 1 paper measuring tape
- 1 direction for use

circaid reduction kit - knee	
Circumference in cm	regular
cE2	< 110
cE	< 102
cC	< 80
one length available:	standard 30cm

Each knee kit includes:

- 1 reduction knee component
- 2 hook and loop knee spines
- 1 paper measuring tape
- 1 direction for use
- 1 reduction kit shelf strap