

SOFT PRO®

HINGED WRIST AIR CONE WHO

PRODUCT FEATURES

- No tools needed
- Provides either positional or progressive extension therapy of the wrist depending on the setting of the hinge
- Effectively opens severe hand and finger contractures with gentle pneumatic stretching
- Lateral forearm support provides optimal control of the forearm, wrist and hand
- Ideal for flexion synergy of the wrist/hand/fingers
- Velcro attached finger separator provides for skin protection at the fingers.
- Color Coded Straps

HCPCS-L3915 (Custom Fitted)
L3916 (Off the Shelf)



INDICATIONS OF USE

The SoftPro Hinged Wrist Air Cone WHO is used to treat flexion synergy of the wrist, wrist/hand pain, contractures of the wrist/hand, Colles fracture, radius fracture, or ulnar styloid fracture positioning