

sealtight® ORIGINAL

ADULT SIZES

ORDER#	SIZE	LENGTH	MAX CIRCUMFERENCE
20100	Hand	12" / 31cm	16" / 41cm
20101	Short Arm	23" / 58cm	16" / 41cm
20102	Long Arm	40" / 101cm	21" / 53.3cm
20103	Short Leg	24" / 61cm	21" / 53.3cm
20104	Long Leg	42" / 107cm	28" / 71cm
20105	Foot / Ankle	12" / 31cm	21" / 53.3cm
20106	Wide Short Leg	24" / 61cm	28" / 71cm
20107	Wide Short Arm	23" / 58cm	21" / 53.3cm

PEDIATRIC SIZES

20200	Small Arm	12" / 31cm	16" / 41cm
20201	Medium Arm	19" / 48cm	16" / 41cm
20202	Large Arm	29" / 74cm	21" / 53.3cm
20203	Small Leg	12" / 31cm	16" / 41cm
20204	Medium Leg	17" / 43cm	21" / 53.3cm
20205	Large Leg	32" / 81cm	21" / 53.3cm

1. Measure cast length.
2. Measure circumference at widest point.
3. Measure circumference around ankle. (heel to heel around ankle bend)

