

Suggested Codes L1686* (Hinge L2624*)

SofTec® Coxa

For post-operative stabilization of the hip.

SofTec Coxa stabilizes the hip joint, for example after total hip replacement or femoral-head resection. It helps prevent dislocation of the hip joint by restricting range of motion with a monocentric hinge.

The hinge is positioned at the height of the anatomical pivot point of the hip joint, and can be adjusted in two planes. It allows for infinite limitation settings of hip abduction/adduction as well as flexion and extension up to 90°. This flexibility helps protect against dislocations and drives positive surgical outcomes.

The orthosis is anatomically shaped and can be easily customized. The pelvic frame, knit parts and tensioning straps can be individually adjusted for optimal fit and function.

The knit fabric around the hip and the thigh support are made of breathable material, and are comfortable on the skin for all day wearing comfort.

Studies show 95.8% of all patients were satisfied with the SofTec Coxa**.

Indications

- Prevention of dislocation
- Total hip-revision surgery
- Arthroscopic FAI procedure
- Other arthroscopic procedures
- Femoral-head resection
- Two and multi-stage revision procedures

*Suggested codes are based on approved coding from PDAC as of 01/21. It is the responsibility of the provider to determine the appropriate billing code, as well as, whether the use of a product complies with medical necessity and other documentation requirements of the payer.

**Observational non-interventional study. Data on file, Bauerfeind



SofTec® Coxa

INNOVATIVE BRACING SOLUTIONS AFTER HIP PROCEDURES



INDIVIDUALIZED CARE,
DOWN TO THE LAST DETAIL.

- STABILIZES AND IMMOBILIZES THE HIP
- RELIEVES PAIN
- SECURES POST-OPERATIVE OUTCOMES



MONOCENTRIC HINGE

- Multidimensional joint: controls range of motion to protect against dislocations
- Extension and flexion can be limited in 10° increments from -10° to +90°
- Adduction and abduction can be adjusted between -6° and +6°



EXCELLENT WEARING COMFORT

- Lightweight, low profile design - can be worn under most pants
- Easy for patients to use and position: can be donned and doffed without help
- High wearing comfort due to soft, breathable fabric
- Customized fit




PELVIC FRAME

- Anatomically shaped pelvic frame for secure positioning of the orthosis on the iliac crests
- Combination of knitted soft fabric and a flexible pelvic shell creates uniform contact pressure to ensure optimal immobilization of the joint



CUSTOMIZATION

- Tool-free adjustment of the hinge
- Height-adjustable thigh belt
- Easily adjustable pelvic/thigh belt and tensioning straps
- Unilateral/bilateral application possible
- Can be adjusted at the bedside on supine patients



Brace: Pelvic + Thigh Support	1	2	3	Order number
Trunk Circumf. in inches	29½–35½	35½–43¼	43¼–51¼	1222521008000°
Trunk Circumf. in cm	75–95	95–115	115–135	

°= size (1-3)

Hinge	Order number
Women (offset)	12225221080001
Men (straight)	12225231080000

For bilateral care - please order an additional hinge and a thigh support.

Thigh Support	Order number
	12925210080000