



INTRODUCING THE GERY...

a comfortable, re-engineered low K2 foot.



The newly re-designed GERY foot aims to provide lower level K2 users with the freedom they deserve.

- Greatly assists your patient's ability to ambulate with it's lightweight design.
- Accommodates heel heights ranging from .2" to .8" (5mm to 20mm).
- Cleaned easily with water and is resistant to dirt thanks to its closed cells EVA foam.



GET BACK TO YOUR LIFE.

- Engineered specifically for the low level K2, AK & BK population. Rated up to 150kg / 330lbs.
- Re-engineered with a new core
- 1 Year warranty
- 100 % Waterproof
- Comfort thanks to its design that reduces residual limb pain. Lightweight thanks to its EVA structure \approx 30% lighter than a SACH foot! Easy step completion thanks to its profiled core.
- Comes with malleolar clip

TECHNICAL FEATURES

| | | |
|--------------------|----------------|---------|
| Order number | 1A200H | 1A200BH |
| Weight | .87lbs. (395g) | |
| Build Height | 3.22" (82mm) | |
| Structure material | EVA | |

| | |
|--------------------------|--------------------------------|
| Heel height (T) | from .2" to .8" (5mm to 20 mm) |
| Connection | Pyramid |
| Maximum patient's weight | 330lbs. (150kg) |
| Activity level | 1-2 |

| REF. | Side | Size | Suggested Code | Color |
|----------|-----------------------|--------------|----------------|-----------|
| 1A200-H | G / D Left / Right | 22cm to 29cm | L5972 | Caucasian |
| 1A200-BH | | | | Ebony |

