



## X2K® Knee Brace



Sizes: s, m, m+, l, xl, xxl

X2K is designed to provide patients with optimal varus and valgus stiffness. The X2K is available in versions with standard or adjustable hinges.

### Common Examples of Use

- ACL
- PCL
- Collateral Ligament protection

### Product Features

- Diamond design for optimal varus and valgus stiffness
- Lightweight Aluminum frame
- Patented adjustable hinges allow for optimal loading across the joint
- Internally mounted straps for ideal brace suspension

Category: ACL and PCL Knee Braces.

### Instructions for Use

[Instructions for Use](#)

### Part Numbers

## Part Numbers

Part #	Description
20001	X2K Custom
20002	X2K Custom w/ Adjustable Hinge
201XX	X2K Prefabricated, Left
202XX	X2K Prefabricated, Right
203XX	X2K Prefabricated, w/ Adjustable Hinge, Left
204XX	X2K Prefabricated, w/ Adjustable Hinge, Right
705XX	X2K Refurbish Kit Kit contains 2 regular condyle pads, 2 thincondyle pads, 2 stop sets, black frame pads, straps, pads, hex key and instructions)

## Product Sizing Chart

X2K / Women's X2K / X2K-OA					
Size #	Size	6" Above Mid Patella		Knee Joint	
XX = 10	XS	13" – 15.5"	(33 – 39 cm)	11" – 12.75"	(30 – 32 cm)
XX = 20	S	15.5" – 18"	(39 – 46 cm)	12.75" – 14.5"	(32 – 37 cm)
XX = 30	M	18" – 19.5"	(46 – 50 cm)	14.5" – 15.25"	(37 – 39 cm)
XX = 35	M+	19.5" – 21"	(50 – 53 cm)	15.25" – 16.25"	(39 – 41 cm)
XX = 40	L	21" – 24"	(53 – 61 cm)	16.25" – 18"	(41 – 46 cm)
XX = 50	XL	24" – 27"	(61 – 69 cm)	18" – 19.75"	(46 – 50 cm)
XX = 60	XXL	27" – 31"	(69 – 79 cm)	19.75" – 22"	(50 – 56 cm)
<i>Circumference should be taken at knee joint and 6" (15 cm) above mid patella</i>					
<i>XS not available in X2K and X2K-OA</i>					