



X2K® OA Knee Brace



Sizes: s, m, m+, l, xl, xxl

Breg X2K-OA knee brace is designed for OA patients who are suffering from osteoarthritis of the knee requiring a brace. The brace provides control of varus / valgus forces associated with unicompartmental osteoarthritis of the knee. Note: X2K OA knee brace requires the use of a hex key (included) for off-loading adjustments.

Common Examples of Use

- Unicompartmental OA-medial
- Chondral defect procedures
- Meniscal procedures

Product Features

- Enhanced adjustable hinge - 12 mm of travel (twice the loading capabilities as the standard adjustable hinge) - both hinges
- Gel comfort condyle pad
- Siliconized strap padding (Straps 1 & 2)
- Over-sized tibial frame pad

Category: Functional OA.

Instructions for Use

[Instructions for Use](#)

Part Numbers

Part Numbers

Part #	Description
251XX	X2K-OA Prefabricated w/ Adjustable Hinge, Left
252XX	X2K-OA Prefabricated w/ Adjustable Hinge, Right

* Refurbish kits include 2 regular condyle pads, 2 gel pads, 2 stop sets, black frame pads with straps, strap tabs, hex key and instructions.

Product Sizing Chart

X2K / Women's X2K / X2K-OA					
Size #	Size	6" Above Mid Patella		Knee Joint	
XX = 10	XS	13" – 15.5"	(33 – 39 cm)	11" – 12.75"	(30 – 32 cm)
XX = 20	S	15.5" – 18"	(39 – 46 cm)	12.75" – 14.5"	(32 – 37 cm)
XX = 30	M	18" – 19.5"	(46 – 50 cm)	14.5" – 15.25"	(37 – 39 cm)
XX = 35	M+	19.5" – 21"	(50 – 53 cm)	15.25" – 16.25"	(39 – 41 cm)
XX = 40	L	21" – 24"	(53 – 61 cm)	16.25" – 18"	(41 – 46 cm)
XX = 50	XL	24" – 27"	(61 – 69 cm)	18" – 19.75"	(46 – 50 cm)
XX = 60	XXL	27" – 31"	(69 – 79 cm)	19.75" – 22"	(50 – 56 cm)

Circumference should be taken at knee joint and 6" (15 cm) above mid patella

XS not available in X2K and X2K-OA