

MODEL 147



Model 147 incorporates the Becker 2041-C hip joint, which allows incremental adjustment of flexion/extension in 20° increments and abduction/adduction in 9° increments.

** All Extra Small models are sent with a variant of our 2040-C EconoHip™ joint. The joints have a preset abduction angle of 15 degrees and are not attached to the pelvic section allowing the Orthotist to choose the position that best suits the needs of the patient.*

Order No.	Size	Pelvic Diameter
147-XS*	X-Small	6" (15cm) to 7-1/2" (19cm)
147-S	Small	7-3/4" (20cm) to 9-1/2" (24cm)
147-M	Medium	9-3/4" (25cm) to 11-1/2" (29cm)
147-L	Large	11-3/4" (30cm) to 13-1/2" (34cm)
147-C	Custom	Custom Central Fab

¹ Banta JV, Scrutton D (2003) *Hip Disorders in Childhood*. Cambridge: Cambridge University Press