

NEW AND IMPROVED TOTAL KNEE® BY ÖSSUR

TOTAL KNEE ENABLES AMPUTEES TO WALK WITH A SMOOTH, MORE NATURAL, ENERGY-EFFICIENT GAIT, ALLOWING FOR INCREASED CONFIDENCE AND STABILITY.

New improvements include:

- Increased durability
- Improved aesthetics
- Easier cosmetic finishing

The Total Knee offers the advantages you've come to rely on:

- Mid-swing shortening provides increased toe clearance to prevent hip-hiking and reduce chances of tripping
- Extension assist and extension promotor ensure full knee extension
- Geometric lock creates stability throughout stance phase
- Adjustable stance flexion bumpers mimic knee flexion during early stance, reducing shock and stress on the residual limb
- Rugged and reliable, yet lightweight





TOTAL KNEE® 2000 | TOTAL KNEE® 1900



Life Without Limitations®

TECHNICAL SPECIFICATIONS



NOTE: Technical specifications are identical for both models



TOTAL KNEE® 2000

For moderately active users who are able, or have potential, to walk at different speeds. 3-phase hydraulic swing control accommodates changes in walking speed. Suitable for moderate activity; walking on level and uneven ground.

- Polycentric knee with geometric locking system
- Swing control adjustments through 3-phase (3-valve) hydraulic system
- Weight limit of 220 lbs (100 kg)
- Adjustable stance flexion
- Adjustable extension promoter

Part Number	Description
2000	Total Knee 2000
2070S	Optional Extension Assist

TOTAL KNEE® 1900

For users with less active lifestyles engaging in low impact activities and walking at a constant speed.

- Polycentric knee with geometric locking system
- Controlled and steady swing phase through use of adjustable friction
- Weight limit of 220 lbs (100 kg)
- Adjustable stance flexion
- Adjustable extension promoter

Part Number	Description
1900	Total Knee 1900
2070S	Optional Extension Assist

SPECIAL KNEE TRIAL PROGRAM: TAKE ONE FOR A TEST DRIVE

To help you make the right choice for your patient, Össur offers a 30-day no obligation trial program on mechanical knees. Speak with your Össur Sales Associate or Customer Service Representative today to request a trial Total Knee 2000 or Total Knee 1900. Once you are sure Total Knee is the right solution for your patient, let us know and we will quickly replace the trial unit with a new Total Knee. Please note that there are no returns of Össur knees once purchased and removed from the box.

ALL TOTAL KNEES HAVE A 24-MONTH WARRANTY.

800 233 6263 FAX 800 663 5982 WWW.OSSUR.COM



Life Without Limitations®