

Measure using patient arm length (elbow to back of the hand).				
Size	Arm Length	Item Number	Optional Left External Rotation Kit	Optional Right External Rotation Kit
S	11.5 - 13 in (29 - 34 cm)	50A10=S	29A222=L-S/M	29A222=R-S/M
M	13 - 14.5 in (35 - 37 cm)	50A10=M	29A222=L-S/M	29A222=R-S/M
L	14.5 - 16 in (38-41 cm)	50A10=L	29A222=L-L	29A222=R-L