

BECKER ORTHOPEDIC SPINAL SYSTEM

MODULAR COMPONENTS FOR HYPEREXTENSION TLSO'S

B.O.S.S. hyperextension TLSO's are completely interchangeable. As the needs of your patient change, different assemblies may be easily added to the existing orthosis. Some orthotists will prefer to stock the modular assemblies, as this method allows great versatility in meeting the needs of the patient. To design your own B.O.S.S. hyperextension TLSO, please follow the steps below. Assemblies may also be ordered individually.

STEP 1: Choose the model number of the thoracic assembly you wish to purchase.

THORACIC ASSEMBLIES

Each set includes two thoracic extensions and sternal pad.



MODEL 25T

Thoracic assembly with swivel sternal pad.

Order No.	Size
25T-XS	X-Small
25T-S	Small
25T-M	Medium
25T-L	Large
25T-XL	X-Large



MODEL 26T

Thoracic assembly with stationary sternal pad.

Order No.	Size
26T-XS	X-Small
26T-S	Small
26T-M	Medium
26T-L	Large
26T-XL	X-Large



MODEL 29T

Thoracic assembly with deltopectoral pads.

Order No.	Size
29T-XS	X-Small
29T-S	Small
29T-M	Medium
29T-L	Large
29T-XL	X-Large

Note: To order Thoracic Assemblies with 15" Boston Extension Uprights, please add the letter "A". For example, to order model 25T-XS with Boston Extensions, your order number would be 25TA-XS.

STEP 2: Choose the model number of the pelvic assembly you wish to purchase.

PELVIC ASSEMBLIES

Each set includes two pelvic extensions and a suprapubic pad.



MODEL 25P

Pelvic assembly with suprapubic pad.

Order No.	Size
25P-XS	X-Small
25P-S	Small
25P-M	Medium
25P-L	Large
25P-XL	X-Large



MODEL 27P

Padded pelvic bar.

Order No.	Size
27P-XS	X-Small
27P-S	Small
27P-M	Medium
27P-L	Large
27P-XL	X-Large



MODEL 28P

Pelvic assembly with anterior and posterior extensions and suprapubic pad.

Order No.	Size
28P-XS	X-Small
28P-S	Small
28P-M	Medium
28P-L	Large
28P-XL	X-Large

STEP 3: Choose the model number of the lateral upright assembly you wish to purchase.


ADJUSTABLE LATERAL UPRIGHTS WITH THORACOLUMBAR PAD

Each set includes two lateral uprights, side pads, thoracolumbar pad, closure system and four lid covers and screws for serrations.

MODEL Q-25AL

Adjustable lateral uprights and thoracolumbar pad with Quick Release closure.


Order No.	Size	Adjustable Range
Q-25AL-S	Small	4-3/4" to 7-1/4"
Q-25AL-M	Medium	6" to 8-1/4"
Q-25AL-L	Large	6-5/8" to 9-3/4"



MODEL L-25AL

Adjustable lateral uprights and thoracolumbar pad with Dot Fastener closure.

Order No.	Size	Adjustable Range
L-25AL-S	Small	4-3/4" to 7-1/4"
L-25AL-M	Medium	6" to 8-1/4"
L-25AL-L	Large	6-5/8" to 9-3/4"



MODEL W-25AL

Adjustable lateral uprights and thoracolumbar pad with Worm Gear closure.

Order No.	Size	Adjustable Range
W-25AL-S	Small	4-3/4" to 7-1/4"
W-25AL-M	Medium	6" to 8-1/4"
W-25AL-L	Large	6-5/8" to 9-3/4"




NON-ADJUSTABLE LATERAL UPRIGHTS WITH THORACOLUMBAR PAD

Each set includes two lateral uprights, side pads, thoracolumbar pad, closure system and four lid covers and screws for serrations.

MODEL L-25L

Non-Adjustable lateral uprights and thoracolumbar pad with Dot Fastener closure.


Order No.	Length (Center of Serration)
L-25L-4.5	4-1/2"
L-25L-5.25	5-1/4"
L-25L-6	6"
L-25L-6.75	6-3/4"
L-25L-7.5	7-1/2"
L-25L-8.25	8-1/4"
L-25L-9	9"
L-25L-9.75	9-3/4"



MODEL W-25L

Non-Adjustable lateral uprights and thoracolumbar pad with Worm Gear closure.

Order No.	Length (Center of Serration)
W-25L-4.5	4-1/2"
W-25L-5.25	5-1/4"
W-25L-6	6"
W-25L-6.75	6-3/4"
W-25L-7.5	7-1/2"
W-25L-8.25	8-1/4"
W-25L-9	9"
W-25L-9.75	9-3/4"



NOTE: Non-Adjustable Lateral Uprights are not available with our Quick Release Closure.


ADDITIONAL MODULAR COMPONENTS



MODEL 25ALE

Adjustable lateral upright assembly.

Order No.	Size	Adjustable Range
25ALE-S	Small	4-3/8" to 7-1/4"
25ALE-M	Medium	6" to 8-1/4"
25ALE-L	Large	6-5/8" to 9-3/4"



MODEL 25ALE-ARM

Adjustable lateral upright slide arm with serration.

Order No.	TLSO Size	Slide Arm Size
25ALE-ARM-XS	XS	Short
25ALE-ARM-S	S	Short
25ALE-ARM-M	M	Short
25ALE-ARM-L	L	Long
25ALE-ARM-XL	XL	Long



MODEL 25-1 STERNAL PAD WITH EXTRA HOLES

Allows greater range of positioning. Please specify small, medium or large when ordering.



MODEL 25-2 SUPRAPUBIC PAD WITH EXTRA HOLES

Allows greater range of positioning. Please specify small, medium or large when ordering.



MODEL 25-3 EXTENDED THORACOLUMBAR PAD

Gives greater support to the thoracolumbar spine. Please specify closure type and small, medium or large when ordering.



MODEL 25-HB BOSTON EXTENSION LATERAL UPRIGHT -15"

ADDITIONAL MODULAR COMPONENTS (not pictured)

Order No.	Description
Q-25-3	Thoracolumbar Pad & Strap for Quick Release Closure (specify size)
Q-25-3A	Dacron Hook & Loop Back Strap for Quick Release Closure (specify size)
Q-25-4A	Latch Arm Buckle & Spring Release Unit for Quick Release Closure
Q-25-5	Metal Base Plate for Quick Release Closure
25-4	Snap-Out Post
25-5	Dot Fastener
25-6	Lid & Lid Screw
25-7	Serrated Disc
25-8	Phillips Head Pad Screw
25-9	Standard Sternal Pad (specify size)

Order No.	Description
25-10	Standard Thoracolumbar Pad (specify size)
25-11	Standard Pubic Pad (specify size)
25-12	Thoracic Extension Arm (specify size)
25-13	Pelvic Extension Arm (specify size)
25-14	Plastic Shield for Lateral Upright (specify size)
25-15	Worm Gear Traction Band - Snap-Out Side
25-16	Worm Gear Traction Band - Fixed Side
25-17	Hook and Loop Strap with DOT Fastener - Long Side (specify size)
25-18	Hook and Loop Strap with Eyelet - Short Side (specify size)
25-19	Plastic Sternal Bar with U-Bracket (specify size)