

## Tuning /ˈtoʊniŋ/ (v) to adjust for optimal performance

LaunchPad designs and develops orthotics on the principles of functionality and comfort. We are on a mission to deliver simple streamlined solutions for customizing motion control and optimizing function. **rotator**, de-rotation straps and **boosters** are two of the newest members of LaunchPads' family of products.

### rotator®

**rotator** de-rotation straps are designed for pediatric patients with internal or external foot rotation that impairs gait. Proper gait mechanics are critical to improving mobility and restoring strength. **rotator**s are made with a latex-free material that combines comfort with the compressive elasticity needed to control rotation. **rotator** can be worn next to the skin or over top of clothing. Suggested L-code L2040 1X

#### Function

- Improve foot projection angle
- Promote heel strike
- Normalize stance phase mechanics
- Promote strengthening

#### Differentiators

- Comfortable compressive material
- No roping
- Used with laces or straps
- Infinite attachment points
- Wraps up from the shoe for easier fitting




### boosters®

**boosters** are a "peel and stack" solution to tune (align) the shank to vertical angle of an AFO to encourage tibial progression over the foot. This time saving product helps to optimize gait performance. **boosters** will never compress and will maintain their position in the shoe thanks to our proprietary 3M adhesive. No more raw material "art projects" in the lab. **boosters** are only meant for AFO's, they are not intended to be used as a heel lift. **boosters** come in 2 Adult sizes and 3 Pediatric sizes. Suggested L-code L3332 1X

Adult sizes (M, L): 1/8" (3.2mm) thick

Pediatric sizes (S, M, L): 3/32" (2.4mm) thick



#### How to Order

See our website [www.LaunchPad-OP.com](http://www.LaunchPad-OP.com) for all our distribution partners.