

LOW ACTIVITY LINER SIZING GUIDELINES

PROXIMAL CIRCUMFERENCE



Measure 15"
from Distal End

DISTAL CIRCUMFERENCE



Measure 1.5"
from Distal End

SIZE	RELAXED	FULLY STRETCHED	RELAXED	FULLY STRETCHED
Medium	9" (23cm)	16" (41cm)	8" (20cm)	12" (30cm)
Medium+	13" (33cm)	21" (53cm)	9" (23cm)	12" (30cm)
Large	15" (38cm)	22" (56cm)	9" (23cm)	15" (38cm)
Large+	16" (41cm)	27" (68cm)	12" (30cm)	16" (41cm)
X-Large	16" (41cm)	30" (76cm)	12" (30cm)	20" (51cm)

All sizes are approximate.