

Locking Hinges

Outside Locking Hinge, Large Heavy Duty

HOSMER

- For an elbow-disarticulation prosthesis
- Five locking positions
- Reciprocating type locking mechanism requires very short cable travel to operate
- Lower anchor strap is laminated around end of humeral socket to insure permanent alignment of hinges
- One locking and one non-locking hinge, with medial side locking

53883 Outside Locking Hinge, *Heavy Duty Large, E-5500R; 15.5oz (434g), Pair, Right*

50950 Outside Locking Hinge, *Heavy Duty Large, E-5500L; 15.5oz (434g), Pair, Left*

53884 Outside Locking Hinge, *Large, No Anchor Strap, E-5500AR; 14.5oz (406g), Pair, Right*

50951 Outside Locking Hinge, *Large, No Anchor Strap, E-5500AL; 14.5oz (406g), Pair, Left*

Includes:

- A. **50966** Locking Bar Assembly
- B. **50820** Spring
- C. **50819** Pin
- D. **50826** Cam
- E. **50822** Keeper
- F. **50821** Driver
- G. **50820** Spring
- H. **50724** Cable Assembly
- I. **50945** Set Screw
- J. **50725** Cable Housing Assembly
- K. **50967** Case
- L. **50958** Long Bushing
- M. **55714** Gear Strap, 50950
53981 Gear Strap, 53883
- N. **53978** Locking Side, 53883
50957 Locking Side, 50950
- O. **50962** Screw
- P. **50959** Screw
- Q. **54140** Free Side Assembly, 53883
50972 Free Side Assembly, 50950
54475 Free Side Assembly, 50951
54476 Free Side Assembly, 53884
- R. **50747** Hanger
- S. **50427** Anchor

