

Locking Hinges

Outside Locking Hinge, Medium

HOSMER


- For an elbow-disarticulation prosthesis
- Seven locking positions
- Reciprocating type locking mechanism requires very short cable travel to operate
- Lower anchor strap is laminated around end of humeral socket to insure permanent alignment of hinges
- Available without lower anchor strap
- Each pair includes one locking and one non-locking hinge, with medial side locking

53879 Outside Locking Hinge, Medium, E-1500R; 7.5oz (210g), Pair, Right

50853 Outside Locking Hinge, Medium, E-1500L; 7.5oz (210g), Pair, Left

53880 Outside Locking Hinge, Medium, No Anchor Strap, E-1500AR; 7.5oz (210g), Pair, Right

50854 Outside Locking Hinge, Medium, No Anchor Strap, E-1500AL; 7.5oz (210g), Pair, Left



Includes:

- A. **50886** Cable Housing Assembly
- B. **50887** Cable Assembly
- C. **50818** Cable Return Spring
- D. **55710** Gear Strap Assembly, 50853 and 50854
- 54607** Gear Strap Assembly, 53879 and 53880
- E. **55583** Locking Bar Assembly
- F. **55709** Case Assembly
- G. **50829** Set Screw
- H. **50884** Cam
- I. **50876** Pin
- J. **50875** Spring
- K. **50879** Keeper
- L. **50878** Driver
- M. **54276** Strap Assembly, 50853
- 50855** Strap Assembly, 53879
- 50856** Strap Assembly, 50854
- 54145** Strap Assembly, 53880
- N. **50808** Screw
- O. **50831** Screw
- P. **50867** Free Side Assembly, 50853
- 54913** Free Side Assembly, 53879
- 50868** Free Side Assembly, 50854
- 54171** Free Side Assembly, 53880
- Q. **50809** Bushing
- R. **50747** Hanger, 50747
- S. **50427** Anchor, 50427

