

Locking Hinges

Outside Locking Hinge, Small

HOSMER

- For an elbow-disarticulation prosthesis
- Seven locking positions
- Reciprocating type locking mechanism requires very short cable travel to operate
- Lower anchor strap is laminated around end of humeral socket to insure permanent alignment of hinges
- Available without lower anchor strap
- Each pair includes one locking and one non-locking hinge, with medial side locking

53881 Outside Locking Hinge, Small, E-2500R; 5oz (140g), Pair, Right


50915 Outside Locking Hinge, Small, E-2500L; 5oz (140g), Pair, Left

53882 Outside Locking Hinge, Small, No Anchor Strap, E-2500AR; 5oz (140g), Pair, Right

50916 Outside Locking Hinge, Small, No Anchor Strap, E-2500AL; 5oz (140g), Pair, Left

Includes:

- A. **50944** Cable Housing Assembly
- B. **50887** Cable Assembly
- C. **55583** Locking Bar Assembly
- D. **54167** Gear Strap Assembly, 53881 and 53882
- 55713** Gear Strap Assembly, 50915 and 50916
- E. **54601** Case Assembly
- F. **50945** Set Screw
- G. **50937** Spring
- H. **50884** Cam
- I. **50876** Pin
- J. **50875** Spring
- K. **50879** Keeper
- L. **50878** Driver
- M. **50920** Strap Assembly, 53881
- 53527** Strap Assembly, 53882
- 50919** Strap Locking Side, 50915
- 53528** Strap Locking Side, 50916
- N. **50933** Bushing
- O. **50932** Screw
- P. **50947** Screw
- Q. **54003** Free Side Assembly, 53881
- 55521** Free Side Assembly, 50915
- 55786** Free Side Assembly, 50916
- 60874** Free Side Assembly, 53882
- R. **50747** Hanger
- S. **50427** Anchor
- T. **51258** Lockwasher



Hinge: 53881