



The Lerman Hip Abduction Rotation Orthosis (HARO) 2000™ from Becker Orthopedic utilizes an anterior pelvic section which provides easier donning and eliminates pressure on the back and spine. Hip joint provides flexion/extension adjustability in 20° increments and abduction/adduction adjustability in 9° increments. The attachment of the joint to the pelvic and thigh modules allows rotational adjustment and accurate alignment in the transverse plane.

The Universal model (145-U) utilizes a low profile pelvic section and a universal thigh section.

Model 145	Small	Medium	Large	X-Large	Plus Size
Waist Circumf	28-32"	34-38"	40-44"	46-50"	52-56"
Left	145-LS	145-LM	145-LL	145-LX	145-LPS
Right	145-RS	145-RM	145-RL	145-RX	145-RPS
Model 145-U			SM/MD	LG/XL	Plus Size
Waist Circumf			28-38"	38-48"	48-58"
Universal			145-U-S-M	145-U-L-X	145-U-PS

Replacement Parts	145	145-U
Full Pelvic Assembly with Liner	145-1	145-U-1
Full Thigh Assembly with Liner	145-2	145-U-2
Pelvic Liner Set	145-5	145-U-5
Pelvic Liner Only	145-6	145-U-6
Posterior Liner/Strap	145-7	145-U-7
Thigh Liner Assembly	145-10	145-U-10