

INTRODUCING

*WillowWood*<sup>®</sup>

LIMBLOGIC<sup>®</sup> 



LimbLogic M is a simple, **elevated vacuum solution** featuring adjustable vacuum and stiffness settings. This sleek, inline, mechanical pump has a low build height and offers clinicians flexibility in componentry selection.

FOR MORE INFORMATION, CALL 800.848.4930  
OR VISIT [WILLOWWOODCO.COM](http://WILLOWWOODCO.COM)

Gary B. - Transtibial Amputee since 2010.



# MECHANICAL ELEVATED VACUUM

## Adjustable Vacuum Settings

For optimal security and patient comfort, clinicians may select Low or High vacuum settings with the turn of an Allen wrench.

## Adjustable Stiffness & Shock Absorption

Adjustable stiffness settings allow LimbLogic M to satisfy a range of patient weights and comfort preferences while reducing impact to a residual limb.

## Standard 4-Hole Mount & Rotatable Distal Pyramid

A 4-hole mount and a rotatable distal pyramid allow for easy, hose-free set-up and alignment.

**BUILD HEIGHT:** 3.6" (9 cm)

**PRODUCT WEIGHT:** 475 g

**U.S. L-CODE SUGGESTIONS:** L5781 (vacuum)  
L5988 (vertical shock)

**TRIAL PERIOD:** 30 days from fit date

**WARRANTY:** 24 months from fit date

### ORDERING INFORMATION:

#### Weight Limits\* by Activity Level

		K2	K3	K4
LLM-1001-L	LimbLogic M, Light	100-170 lb (45-77 kg)		
LLM-1000-L	LimbLogic M, Standard	160-330 lb (73-150 kg)	160-300 lb (73-136 kg)	

*\* Body weight plus any loads normally or routinely carried*

### Related Products:

LLM-KIT-L	LimbLogic M Fabrication Kit	LLM-11001	Diagnostic Kit
-----------	-----------------------------	-----------	----------------

Refer to [willowwoodco.com](http://willowwoodco.com) for video resources and instructions for use.

