

Socket Comfort – Silicone Liners



width of textile strip: 38 mm

Dimension sketch:



total length: approx. 370 mm



AK-Control.Sil

Trans-femoral silicone liner made of medium stable Control-Sil-Silicone, with distal connection and anti-rotation strips

Designated use

- activity level (walky): 1 to 4
- ideal for liner-care of above-knee amputations and for care of voluminous below-knee amputations

Characteristics

- 2 mm wall thickness for controlled expansion in the socket entry zone
- outer matrix (sewed-on strip) for reduction of longitudinal expansion (inhibited longitudinal expansion)
- anti-rotation strips on the outer matrix prevent twisting between liner and prosthetic socket
- especially soft functional textile cover for easily rolling the liner up or down
- very good radial expansion adjusts to volume changes of the residual limb
- rotation control to improve guidance of prosthesis
- with distal connection for all Streifeneder pin types and locking systems

Material

- wall thickness: 2 mm, see sketch
- Control-Sil-Silicone, medium stable, antibacterially fitted
- approx. 40 Shore 00
- abrasion-proof textile cover made of polyamide/elastane
- three rows of silicone stripes on the textile strip in vertical direction
- with PUR-rotation control
- colour: silver grey / black

Warranty

- 6 months

Size	Wall thickness	Item-No.
25, 26,5, 28, 30, 32, 34, 36, 38, 40, 42, 45, 50, 55	2 mm	60L1/..

Ordering example: 60L1/25 = size 25. Scope of supply does include stand.



Sizing: Measure the distal circumference approx. 4 cm above the stump end and select the liner in the same size or one size smaller.