

ANTERIOR POSTERIOR TAPERED

APDT-TP | APFR-TP


Gel Liner —
With or Without Distal Attachment



DESCRIPTION

The ALPS™ AP Tapered Liner is gradually tapered from the anterior to the posterior to provide superior comfort. This liner is available in a pin and lock system or suction suspension. The AP Tapered Liner features our black high performance fabric with gel to assist in reducing bunching in the popliteal region.

TOP FEATURES




► Constructed using high performance black fabric with gel to reduce bunching in the popliteal region

► Tapered - 6mm Anterior to 3mm Posterior

CROSS SECTION DIAGRAM


6.0 MM TO 3.0 MM TAPERED



HOW TO SIZE

Measurements taken
6cm from the distal end

SIZING CHART



MEASURED SIZE (CM)	LOCKING 6.0 MM or 3.0 MM UNIFORM	CUSHION 6.0 MM or 3.0 MM UNIFORM
20-23 cm	APDT 20 - TP	APFR 20 - TP
24-25 cm	APDT 24 - TP	APFR 24 - TP
26-27 cm	APDT 26 - TP	APFR 26 - TP
28-31 cm	APDT 28 - TP	APFR 28 - TP
32-37 cm	APDT 32 - TP	APFR 32 - TP



T: 727.528.8566 | 800.574.5426 | F: 727.528.8862

www.easyliner.com | info@easyliner.com

2895 42nd Avenue North St. Petersburg, Florida 33714

