



SPL2

from Fillauer LLC

The most popular
stance control
orthotic knee joint
just got better



- Enhanced position dependent pendulum design
- Multiple locking positions for tentative gait
- Added free motion feature
- Easy fabrication with minimal follow-up maintenance
- Works with any orthotic ankle design
- Only design available for use as a rehab knee orthosis
- Bilateral or unilateral uprights for easy donning

*Available in a kit or as a Fillauer
Custom Fabrication orthosis*

U.S. Patent
6979304



Fillauer LLC

fillauer.com/products/spl2

1.800.251.6398

AD278/07-09/08-18-10

Swing Phase Lock 2 Orthotic Knee Joint Ordering Guide

INDICATIONS

- Post-polio
- Spinal involvement
- Cerebral vascular accident
- Peripheral paresis/paralysis
- Nerve inflammations
- Neurological failures
- Myopathies
- Multiple sclerosis (or similar diseases)

CONTRAINDICATIONS

- Knee flexion contractures > 10°
- Central paresis/paralysis
- Hip flexion contracture
- Weight bearing orthotic applications

WARRANTY

Fillauer orthotic components are warranted for 12 months from the date of shipment from the Fillauer warehouse or from an authorized Fillauer distributor. Items under warranty will be replaced or repaired (at Fillauer's discretion) at no charge. The warranty will be void if the item has been fabricated or installed outside Fillauer's recommendations, if the item has been exposed to a corrosive environment, or if the item has been used in extremely abusive activities that could result in injury.

SWING PHASE LOCK KITS

IQ Swing Phase Lock 2 Kit includes the SPL2 Joint, SPC Joint, and Satellite (remote control)

IQ200L	IQ Swing Phase Lock 2 Kit, Left
IQ250L	IQ SPL2 Joint, Left
IQ120L	IQ SPC Joint, Left
IQ140	IQ SPL Satellite (remote control)
IQ050L	IQ Single Sided Swing Phase Lock 2
IQ250L	IQ SPL2 Joint, Left
IQ140	IQ SPL Satellite (remote control)
IQ200R	IQ Swing Phase Lock 2 Kit, Right
IQ250R	IQ SPL2 Joint, Right
IQ120R	IQ SPC Joint, Right
IQ140	IQ SPL Satellite (remote control)
IQ050R	IQ Single Sided Swing Phase Lock 2
IQ250R	IQ SPL2 Joint, Right
IQ140	IQ SPL Satellite (remote control)


NEW FEATURES OF THE SPL2

Two Locking Positions at 0° and 15°

PEEK Bushing Insert

Extended Wear Cylindrical Bumper

Teflon Bushings for Sideloads



SWING PHASE LOCK KITS

Fabrication Kit

IQ150 Alignment Kit

Upright Kits

023516 Swing Lock Upright Kit, AL 3/16 x 3/4

023616 Swing Lock Upright Kit, SS 3/16 x 3/4

Polycar-C™

076630 PolyCar-C™ Insert, .060 Small, pair

076632 PolyCar-C™ Insert, .090 Med, pair

076634 PolyCar-C™ Insert, .090 Large, pair

076636 PolyCar-C™ Insert, .125 XL, pair

020803 Pre-cut PolyCar-C™, Large

Carbon also available

Ankle Joints

020886 Double Action, SS, Large, 3/4 channel

020888 Double Action w/Positioning Pins, SS, Large, 3/4 channel

020882 Double Action, SS, Large, 5/8 channel

020884 Double Action w/Positioning Pins, SS, Large, 5/8 channel

Insert Stirrups for Plastic Orthotics

020099 Insert Stirrup, Limited Motion, 2-5/16"

020115 Insert Stirrup, Limited Motion, 2-3/4"

020123 Insert Stirrup, Toe Lift, 2-5/16"

020149 Insert Stirrup, Toe Lift, 2-3/4"

020094 Insert Stirrup, Double Action, 2-3/4"

020016 Offset Stirrup, Limited, 2-1/4" x 3/8"

020024 Offset Stirrup, Limited, 2-1/4" x 1/2"

020032 Offset Stirrup, Limited, 2-1/2" x 3/8"

020040 Offset Stirrup, Limited, 2-1/2" x 1/2"

020057 Offset Stirrup, Toe Lift, 2-1/4" x 3/8"

020065 Offset Stirrup, Toe Lift, 2-1/4" x 1/2"

020073 Offset Stirrup, Toe Lift, 2-1/2" x 3/8"

020081 Offset Stirrup, Toe Lift, 2-1/2" x 1/2"